VOCO

Bifix® Veneer LC

Instructions for use

MD EU Medical device

In accordance with DIN EN ISO 4049

Product description:

Bifix Veneer LC is a light-curing, radiopaque nanohybrid luting composite, resistant to discoloration, for the adhesive luting of translucent veneers, inlays and onlays made of ceramic or composite. **Bifix Veneer LC** contains 81% by weight inorganic fillers in a light-curing methacrylate matrix. The specially formulated thixotropic properties of **Bifix Veneer LC** allow simple placement of the restoration. The product has excellent flow properties with simultaneous setting of the excess material, which is therefore easy to remove, e.g. after brief polymerisation (tack curing).

Shades:

Bright, Clear, Warm

Indications:

Permanent luting of ceramic and composite veneers, inlays and onlays in a layer thickness of not more than 2 mm with sufficient translucency of the restoration.

Contraindications:

Bifix Veneer LC contains methacrylates and BHT. **Bifix Veneer LC** should therefore not be used in patients with a known hypersensitivity (allergy) to these ingredients.

Patient target group:

Bifix Veneer LC is suitable for use in all patients without any age or gender restrictions.

Performance features:

The product's performance features satisfy the requirements of the intended use and the relevant product standards.

User:

Bifix Veneer LC should only be used by a professionally trained dental practitioner.

Use:

Bifix Veneer LC must be applied using the adhesive technique. If stored in the fridge, the materials must be brought up to room temperature before use. The **Sticky Tim** application aid can be used to improve handling of the restoration.

1. Shade selection (optional):

- 1.1. Remove the temporary restoration, if present. Clean the tooth, rinse thoroughly and dry gently with air. Avoid overdryng.
- 1.2. The shade appearance of the cured **Bifix Veneer LC** can be simulated by placing the permanent restoration with **Bifix Veneer Try-In** in the intended shade. To do this, lock the supplied type 50 application tip on the syringe by turning it clockwise and apply **Bifix Veneer Try-In** to the adhesive surface of the restoration before carefully positioning it on the tooth.
- 1.3. Following the shade check, remove the restoration and rinse the **Bifix Veneer Try-In** off the tooth and restoration thoroughly with a water spray, and then dry both. Take care not to overdry the dentine. Any residual **Bifix Veneer Try-In** can affect the permanent luting of the restoration.

2. Pre-treatment and conditioning of the restoration:

For detailed information on correct use, also follow the instructions provided with the bonding agent used as well as the manufacturer's instructions for preparing the particular restorations.

2.1. Preparation of the restoration

- a) For the preparation of silicate ceramic, the adhesive surface of the restoration being placed must be conditioned with a hydrofluoric acid etching agent. To do this, prepare it with a suitable acid-etch agent (following the manufacturer's instructions) and then rinse thoroughly with water before drying with air.
- b) For the preparation of zirconium dioxide or composite: If not already done at the dental lab, sandblast the adhesive surface with a suitable abrasive material (follow the manufacturer's instructions). Carefully remove any residual abrasive material, e.g. using an ultrasound bath, and then dry with air. Do not use phosphoric acid.

- 2.2. Apply **Ceramic Bond** (or other suitable bonding agent; follow the instructions for use) to the adhesive surface of the prepared restoration and allow it to work for 60 s. Then dry the surface carefully with air. Do not touch or otherwise contaminate the conditioned surface.

3. Isolation and moisture control:

Ensure that the work area is sufficiently dry. Furthermore, interdental matrix strips should be used to avoid undesirable bonding to the neighbouring teeth.

4. Adhesive:

- The instructions provided by the adhesive manufacturer must be followed. These instructions describe the process for using the adhesive **Futurabond U**.
- 4.1. Clean the cavity or preparation site, for example with pumice stone paste, then rinse and dry off excess water with a gentle air stream. Do not overdry the dentine. The objective is a surface with a silk matt appearance.
 - 4.2. Etching of the dental hard tissue (recommended)
Futurabond U ensures a secure bond to the natural dental hard tissue in both self-etch mode and total etch mode. Etching can additionally boost the adhesive strength on enamel especially. Unprepared enamel must always be etched at the beginning of the bonding process. Etch the surfaces using a dental acid-etch agent based on phosphoric acid, e.g. **Vococid** (see applicable instructions for use). Aspirate the acid-etch agent and rinse with water. Dry off excess moisture with a gentle stream of air. Do not overdry the dentine. The objective is a surface with a silk matt appearance.
 - 4.3. Application of **Futurabond U**

Apply the adhesive evenly to the tooth and rub in for 20 s using the **Single Tim**. Dry the adhesive layer with dry air for at least 5 s in order to remove the solvent. Cure using a standard light-curing unit for 10 s (LED or halogen lamp with a light output of at least 500 mW/cm²).

5. Placing the restoration with Bifix Veneer LC:

- 5.1. Lock the supplied type 50 application tip on the syringe by turning it clockwise and apply **Bifix Veneer LC** directly. Apply a sufficient layer of the selected shade of **Bifix Veneer LC** to the conditioned adhesive surface of the prepared restoration using the application tip. Avoid direct exposure of the material to the treatment light. To make sure that the syringe does not drip or run, the plunger must not be retracted during or after use. The attached application tip is intended for single use only. Remove it after use and seal the syringe again tightly with its original cap. Protect the syringe from contamination.
- 5.2. Place the restoration quickly, exerting moderate pressure so that excess material escapes on all sides.
- 5.3. Fix the restoration in place while removing all the excess material (e.g. by light-curing a spot in the centre of the restoration). Select one of the following techniques for removal of the excess material:

- a) **Wiping technique:** Remove excess material with a disposable brush, foam pellet and floss.
 - b) **Tack-Curing:** Cure excess material along the edges of the restoration very briefly with a conventional light-curing unit held at a distance of a few centimetres so as to achieve a gel-like state. The excess is then easy to remove using a probe or scaler. A good consistency for removal of the excess material depends on the light-curing unit used. Take special care to remove excess material in hard-to-reach areas completely.
6. **Polymerisation:**
Conventional curing lights with a wavelength in the range of 400–500 nm are suitable for the final light-curing of the material. Light-cure the restoration one segment at a time from every accessible surface, beginning with the approximal areas.
The minimum curing time is 20 s at 500–1,000 mW/cm² or 10 s at $\geq 1,000$ mW/cm².
In order to minimise inhibition by atmospheric oxygen around the edges of the restoration, they should be covered with a transparent glycerine gel prior to final curing. Once curing is complete, the glycerine gel can be rinsed off with water.
7. **Finishing:**
Remove the matrix strips and carefully finish and polish the edges of the restoration with instruments suitable for composite. Check the occlusion.

Warnings, precautionary measures:

- Phenolic substances, especially products containing eugenol and thymol, interfere with the curing of filling composites. For this reason, the use of zinc oxide eugenol cements or other materials containing eugenol in combination with **Bifix Veneer LC** should be avoided.
- If adhesives or silane bonding agents other than **Ceramic Bond** and **Futurabond U** are used with **Bifix Veneer LC**, the instructions for use of the relevant products must be followed.
- Our information and/or advice do not relieve you of the obligation of checking that the products supplied by us are suitable for their intended purpose.

Constituents (in descending order according to content):

Barium aluminium borosilicate glass, silicon dioxide, HEDMA, BisGMA, TEGDMA, BisEMA, fumed silica, initiators, stabilisers, pigments

Storage instructions and application method:

Storage at 4°C–28°C. Seal the syringes immediately after dispensing the material to avoid exposure to light and resultant polymerisation. Do not use after the expiry date.

Disposal:

Dispose of the product in accordance with local regulations.

Reporting obligation:

Serious events such as death, temporary or permanent serious deterioration of a patient's, user's or other person's health condition and a serious risk to public health that arise or could have arisen in association with the use of **Bifix Veneer LC** must be reported to VOCO GmbH and the responsible authority.

Note:

The Summary of Safety and Clinical Performance of **Bifix Veneer LC** can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>). Detailed information can also be found at www.voco.dental.

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