



VPS impression material

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
PUTTY SOFT / PUTTY HARD - Fast Set
TRAY MATERIAL / TRAY MATERIAL MAXI

LIGHT BODY - Fast Set
ULTRA LIGHT - Fast Set

Proclinic Expert	Light Body Regular Set	Regular Body Regular Set	Light Body / Ultra Light Body Fast Set	Heavy Body Regular Set	Heavy Body Fast Set	Putty Soft Regular Set	Putty Soft / 380 Putty Soft Fast Set	Putty Hard Fast Set
Reorder	15153	15152	15154 / 55/ 56	15158	15159	15150	15151 / 15157	15163
ISO 4823 ADA 19	Type 3 Light-bodied Consistency	Type 2 Medium- bodied Consistency	Type 3 Light-bodied Consistency	Type 1 Heavy- bodied Consistency	Type 1 Heavy- bodied Consistency	Type 0 Putty Consistency	Type 0 Putty Consistency	Type 0 Putty Consistency
Dosierung Base: Katalysator	1:1	1:1	1:1	5:1	5:1	1:1	1:1 / 5:1	1:1
Anmischzeit	Auto mixing System	Auto mixing System	Auto mixing System	Auto mixing System	Auto mixing System	30" Auto mixing System	30" Auto mixing System	30" Auto mixing System
Gesamtverarbeitungszeit*	2'	2'	1'30"	2'	1'30"	2	1'30"	1'30"
Verweildauer im Mund (min.)	3'30"	3'30"	2'30"	3'30"	2'30"	3'30"	2'30"	2'30"
Abbindezeit*	5'30"	5'30"	4'	5'30"	4'	5'30"	4'	4'
Verformung unter Druck (Min - Max)	3 - 5%	3 - 5%	3 - 5%	3 - 5%	3 - 5%	1 - 3%	1 - 3%	1 - 3%
Rückstellung nach Verformung	> 99,5%	> 99,5%	> 99,5%	> 99,5%	> 99,5%	> 99,0%	> 99,0%	> 99,0%
Lineare Dimensionsabweichung nach 24 Stunden	< 0,20%	< 0,20%	< 0,20%	< 0,20%	< 0,20%	< 0,20%	< 0,20%	< 0,20%
Härte "Shore A"	45±2	45±2	45±2	60±2	60±2	60±2	60±2	70±2

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HYDROPHILIC ADDITION SILICONE (VINYLPOLYSILOXANE) SUITABLE FOR HIGH PRECISION IMPRESSIONS

USE

Proclinc Putty Soft: Vinylpolysiloxane for first impressions with elastic final consistency, recommended for all impression techniques: two-stage(double impression*) or single-stage (simultaneous).

Proclinc Putty Hard: Vinylpolysiloxane for first impressions with high final hardness, recommended for two-stage(double impression*) techniques.

Proclinc Tray Material: High viscosity vinylpolysiloxane. Recommended for filling impression trays in conjunction with a 'wash' within the Proclinc range (i.e. light body, extra light body & regular body) in the single-stage (simultaneous) technique.

Proclinc Regular Body: Medium viscosity hydrophilic vinylpolysiloxane. Recommended for use in conjunction with Proclinc Putty or Tray Material in the single-stage (simultaneous) technique, or as a monophasic material in a special tray.

Proclinc Light Body: Low viscosity hydrophilic vinylpolysiloxane. Recommended for use in conjunction with Proclinc Putty or Tray Material in the two-stage (double impression*) or single-stage (simultaneous) technique.

Proclinc Ultra Light Body: Very low viscosity hydrophilic vinylpolysiloxane. Recommended for use in conjunction with Proclinc Putty or Tray Material in the two-stage (double impression*) or in the single-stage (simultaneous) technique.

*When using this product in the two-stage (double impression) technique, we recommend that you place a polythene spacer sheet on the putty during the first stage of impression taking to improve detail quality in the next stage.

IMPRESSION TRAYS: PREPARATION AND ADHESIVES

All standard and/or special trays can be used, provided they are rigid and non-deformable. The range of **Proclinc's** impression trays are especially recommended.

Proper adhesion of the material to the impression tray is essential to avoid errors, especially during removal from the mouth. To ensure a strong bond, make sure that the adhesive is appropriate for the type and material of the tray.

IMPORTANT: when using special trays, ask the technician if the resin used is self or light curing.

For correct use, follow the instructions for each individual adhesive.

Proclinc PUTTY

MIXING:

- Take the 2 components out of the tubs using the colour coded dosing spoons.
- Mix identical proportions of Base and Catalyst; until you obtain a homogeneous mix, with no streaks (about 30 sec.). Then place into the impression tray.

WARNINGS:

- Proportions must be strictly observed in order to have the right setting times: a larger quantity of catalyst does not accelerate the setting time.

- Close tubs immediately after use. Do not invert lids, lower lids and dosing spoons.
- Contact with latex gloves may prevent polysiloxanes from hardening. Avoid direct or indirect contact with latex gloves. Instead, wash your hands well to eliminate all traces of impurities or use vinyl gloves.
- Take care not to stain your clothing as vinyl polysiloxanes are chemically resistant.

Proclinic FLUIDS CARTRIDGE VERSION

PREPARING THE DEVICE:

Assemble the device as shown in the drawings (1-2).

Warning: Before assembling the tip, ensure that the two components (base and catalyst) flow out evenly by applying a light pressure on the dispenser lever and extruding a small amount of material, which must be removed (3).

Subsequently, insert the mixing tip and, if required, the intraoral tip in the cartridge. (4)

Proclinic Putty, Tray Material (Maxi Version 380ml, 5:1)

PREPARING THE DEVICE

- Press in the levers and remove the protective cap (fig. 5).
- Insert the cartridge into the machine as instructed by the manufacturer.
- Ensure an equal length of material is coming out. Clean and remove any excess material (fig. 6).
- Fit and secure the dynamic-static mixing tip (fig. 7).
- Fit the **mixing tip locker** taking care to engage the clips (fig. 8A) on the cartridge (fig. 8B).
- The Proclinic 5:1 system cartridge is now ready to use (fig. 10).
- Leave the dynamic-static mixing tip in place after use.
To replace it, simply turn the **mixing tip locker** through 45° and remove it (fig. 11A). Then remove the dynamic-static mixing tip (fig. 11B).
- Follow the machine manufacturer's instructions to remove and replace the cartridge.

WARNINGS

- Always dispense equal amounts of base & catalyst before fitting the dynamic-static mixing tip. If you forget to do so, start the machine and wait until the material coming out is uniform in colour, then take away any non-homogenous material before proceeding.
- Only use original Proclinic 5:1 system cartridges, dynamic-static mixing tips and **mixing tip lockers**. The use of other mixing tips does not guarantee optimal mixing and can jeopardise the product's characteristics and end results.
- The **mixing tip locker** cannot be used with cartridges not from Proclinic's 5:1 mixing system.

WARNING: In sensitive subjects, vinyl polysiloxanes may cause irritation or other allergic reactions.

CLEANING AND DISINFECTION OF THE IMPRESSION

After rinsing well under running water, the impression can be immediately disinfected.

CASTING MODELS

There is no need to wait; models may be cast immediately after disinfection, or up to two weeks later, provided that the impression is kept at room temperature. Proclinic is compatible with all the best plasters on the market and with polyurethane resins. We recommend Proclinic gypsum, type 4 die stone.

CLEANING AND DISINFECTING THE DISPENSER

The dispenser can be cleaned, disinfected, and sterilised by immersion in liquid or in a steam autoclave.

STORAGE

Impressions made with Proclinic must be kept at room temperature. Proclinic impression materials are guaranteed for 3 years following the date of production provided they are kept at a temperature of between 5° and 27° C / 41°- 80°F.

GALVANISATION

Proclinic impressions may be galvanised with copper or silver.


IMPORTANT REMARKS: Consulting services on the use of our products, provided in any form, be that verbal, written or through demonstrations, are based upon the current state of dentistry and of our know-how. They are to be considered as non-binding information, also in relation to possible third party claims, and do not exonerate the user from personally checking whether the product is suitable for the planned use. Use and application by the user take place under no possible control by the Company, and are therefore under the user's responsibility. Possible liability for damages is limited to the value of goods supplied by the Company and utilized by the user.

FOR DENTAL USE ONLY

TIMES FOR CLINICAL USE

PUTTY SOFT / PUTTY HARD - Fast Set
TRAY MATERIAL / TRAY MATERIAL MAXI

LIGHT BODY - Fast Set
ULTRA LIGHT - Fast Set

Proclinic Expert	Light Body Regular Set	Regular Body Regular Set	Light Body / Ultra Light Body Fast Set	Heavy Body Regular Set	Heavy Body Fast Set	Putty Soft Regular Set	Putty Soft / 380 Putty Soft Fast Set	Putty Hard Fast Set
Reorder	15153	15152	15154 / 55/ 56	15158	15159	15150	15151 / 15157	15163
ISO 4823 ADA 19	Type 3 Light-bodied Consistency	Type 2 Medium- bodied Consistency	Type 3 Light-bodied Consistency	Type 1 Heavy- bodied Consistency	Type 1 Heavy- bodied Consistency	Type 0 Putty Consistency	Type 0 Putty Consistency	Type 0 Putty Consistency
Mixing ratio Base:Catalyst	1:1	1:1	1:1	5:1	5:1	1:1	1:1 / 5:1	1:1
Mixing time	Auto mixing System	Auto mixing System	Auto mixing System	Auto mixing System	Auto mixing System	30" Auto mixing System	30" Auto mixing System	30" Auto mixing System
Total working time*	2'	2'	1'30"	2'	1'30"	2	1'30"	1'30"
Time in mouth (minimum)	3'30"	3'30"	2'30"	3'30"	2'30"	3'30"	2'30"	2'30"
Setting time*	5'30"	5'30"	4'	5'30"	4'	5'30"	4'	4'
Strain in compression (min - max)	3 - 5%	3 - 5%	3 - 5%	3 - 5%	3 - 5%	1 - 3%	1 - 3%	1 - 3%
Recovery from deformation	> 99,5%	> 99,5%	> 99,5%	> 99,5%	> 99,5%	> 99,0%	> 99,0%	> 99,0%
Linear dimensional variation after 24h	< 0,20%	< 0,20%	< 0,20%	< 0,20%	< 0,20%	< 0,20%	< 0,20%	< 0,20%
Hardness "Shore A"	45±2	45±2	45±2	60±2	60±2	60±2	60±2	70±2

*The times mentioned must be intended from the start of mixing at 23°C - 73°F. Higher temperatures reduce the times, lower temperatures increase them.

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
LIGHT BODY - Fast Set

ULTRA LIGHT - Fast Set

Proclinic Expert	Light Body Regular Set	Regular Body Regular Set	Light Body / Ultra Light Body Fast Set	Heavy Body Regular Set	Heavy Body Fast Set	Putty Soft Regular Set	Putty Soft / 380 Putty Soft Fast Set	Putty Hard Fast Set
Reorder	15153	15152	15154 / 55/ 56	15158	15159	15150	15151 / 15157	15163
ISO 4823 ADA 19	Type 3 Light-bodied Consistency	Type 2 Medium- bodied Consistency	Type 3 Light-bodied Consistency	Type 1 Heavy- bodied Consistency	Type 1 Heavy- bodied Consistency	Type 0 Putty Consistency	Type 0 Putty Consistency	Type 0 Putty Consistency
Dosage Base: Catalyseur	1:1	1:1	1:1	5:1	5:1	1:1	1:1 / 5:1	1:1
Temps de mélange	Auto mixing System	Auto mixing System	Auto mixing System	Auto mixing System	Auto mixing System	30" Auto mixing System	30" Auto mixing System	30" Auto mixing System
Temps total de travail*	2'	2'	1'30"	2'	1'30"	2	1'30"	1'30"
Temps de prise en bouche (min)	3'30"	3'30"	2'30"	3'30"	2'30"	3'30"	2'30"	2'30"
Temps de prise*	5'30"	5'30"	4'	5'30"	4'	5'30"	4'	4'
Déformation sous compression	3-5%	3-5%	3-5%	3-5%	3-5%	1-3%	1-3%	1-3%
Mémoire élastique	> 99,5%	> 99,5%	> 99,5%	> 99,5%	> 99,5%	> 99,0%	> 99,0%	> 99,0%
Variation dimensionnelle linéaire après 24 heures	< 0,20%	< 0,20%	< 0,20%	< 0,20%	< 0,20%	< 0,20%	< 0,20%	< 0,20%
Dureté "Shore A"	45±2	45±2	45±2	60±2	60±2	60±2	60±2	70±2

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