bestdent

CE
0197

ES

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Content

1 Product introduction.....	38
2 Operation instructions.....	46
3 Cleaning, disinfection and sterilization	49
4 Maintenance.....	64
5 Storage and transport conditions	64
6. Environment protection.....	65
7 After-sales service	66
8 Symble.....	66
9 Statement	69

1 Product introduction

1.1 Overview

It is used to clamp rotating instruments such as burs, dental drills, and dental files, and drive their movements to realize dental operations such as cutting, grinding and drilling. Features:

- a) Low noise.
- b) Adjustable accurate speed; stable and strong cutting force.
- c) Long service life.
- d) Adaptable to electric motor that has standard E shape connector.

1.2 Technical specifications

Model	WJ-11L	WJ-11
Chuck Type	Gland chuck	

Gear ratio	1:1
Connector type	Standard E shape interface ((ISO 3964-2016))
Type of rotating equipment	ISO 1797-1-2011,type 1,φ2.35mm needle
Maximum length of the rotary instrument	22.5mm
Maximum working diameter of the rotary	4mm
Rod type and size	ISO 1797-1-2011, type1 φ2.35mm
Minimum suitable length of the rod	12.5mm
Spray water pressure setting range	0.5-2 (recommended 2) bar
Spray air pressure setting range	1.5-3 (recommended 2) bar

Model	WJ-11L	WJ-11
Spray water consumption at 2 bar water pressure	> 50ml/min	
Spray gas consumption at 2 bar gas pressure	> 1.5NI/min	
Lighting method	Fiber optic	-
Recommended light source	LED lamp (light level: 20000-65000lx)	-

Note: WJ-11L with optical fiber can guide light, WJ-11 can not guide light, in addition to the difference, there is no difference in the use and maintenance.

1.3 Intended use

It is used to clamp rotating instruments such as burs, dental drills, and dental files, and drive their movements to realize dental operations such as cutting, grinding and drilling.

1.3.1 Intended user: Professional dentists.

1.3.2 Intended patient Patients who need dental treatment procedure for the removal of decayed matter, cavity and crown preparations, and removal of fillings and surface finishing of tooth and restoration surfaces

1.4 Structure and compositions

Contra-angle Handpiece mainly consists of contra-angle, driving system, handle, etc. these components can not site disassembly.

1.5 Configuration

Please refer to packing list for more details of configuration.

1.6 Operation conditions

1.6.1 Environment temperature: +5° C~+40° C

1.6.2 Relative humidity: 30%~75%

1.6.3 Atmospheric pressure: 70kPa~106kPa

1.6.4 Cooling water temperature: $\leq +25^{\circ}$ C

1.7 Device safety classification

1.7.1 Type of protection against electric shock: class II equipment (when conneted to the Dental Electric Motor).

1.7.2 Degree of protection against electric shock: B type applied part (when conneted to the Dental Electric Motor).

1.7.3 Classified by operation mode: Intermittent operating device (when conneted to the Dental Electric Motor).

1.7.4 Degree of safety application in the presence of a flammable anesthetic mixture with air, oxygen, or nitrous oxide: Equipment cannot be used in the presence of a flammable anesthetic mixture with air, oxygen, or nitrous oxide. 2.5.6 The equipment cannot be used in the nuclear magnetic resonance environment.

1.7.5 The equipment cannot be used in the nuclear magnetic resonance environment.

1.8 Warnings

1.8.1 Please check the operation outside oral cavity before treatment. Pay full attention to check whether there is vibration, noise, or heating. If there is any unexpected problem, please stop operation.

1.8.2 Avoid strong impact and dropping;

1.8.3 Do not arbitrarily disassemble or modify the handpiece;

1.8.4 Please prepare a backup Contra-angle Handpieces to cope with the case that handpiece malfunction during operation;

- 1.8.5 Do not use non-recommended needles. Otherwise, the running needle may bend or flyout;
- 1.8.6 Please keep the needle handle clean. Otherwise, the stains on the surface of handle will affect the correct fixation of needle, causing the needle to fly out during operation;
- 1.8.7 Cooling water should be supplied during operation. Otherwise it may cause overheating;
- 1.8.8 For safety, please wear protective goggles or dust mask while cutting;
- 1.8.9 Please carry out regular maintenance and inspection;
- 1.8.10 If the contra-angle has not been used for a period of time, before operation, please maintain it as per Clause 4 Maintenance.
- 1.8.11 Damage or serious wear will lead to increased transmission resistance, abnormal sound and so on. As a result, the vibration of Contra-angle Handpieces will be strong; the drill will shake greatly; and there will be abnormal sound, etc. The operator should timely replace the contra-angle according to clinical conditions. If the contra- angle wears to a certain extent, please replace that with a new one.

1.8.12 Improper use will damage this product, and then hurt the patients, users and the third parties.

1.8.13 Any serious incident that has occurred in relation to the device should be reported to manufacturer and the competent authority of the Member State in which the user and/or patient is established.

1.9 Risk information

1.9.1 Clinical benefits Unknow.

1.9.2 Contraindications Unknow.

1.9.3 Residual risks and undesirable side-effects Unknow.

1.10 Safety precautions

1.10.1 For repairs and purchase of spare parts, please contact our authorized supplier.

1.10.2 The accuracy of the speed monitoring depends on the high-precision performance of the handpiece installed on the micro motor. If the handpiece of other manufacturers is used, the actual speed value may not be displayed correctly. To ensure the actual matching display speed, please use the

matching dental electric motor.

1.10.3 Read this operating manual Handpiece before use and fully understand the functions of each part and connect method.

1.10.4 Check the Dental Low Speed before use to confirm that there is no abnormality. If the device is permanently malfunctioning (excessive vibration, noise and heat generation, etc.), please immediately close it and return it to the authorized dealer.

2 Operation instructions

2.1 Installing/Removing the handpiece

2.1.1 Installation Insert the motor plug into the contra angle, when connecting with the fiber optic contra angle, align the locking nitch on the contra angle with the positioning notch on the motor, if the contra angle is installed in place, a "click" can be heard.

2.1.2 At this time there should be no relative rotation between the contra-angle handpiece and the motor.

If no “click” is heard, rotate the contra angle until the sound is heard, at which point the contra angle is not rotated relative to the motor.

Disassemble the motor and contra angle by gripping them and pulling them out directly, as in Figure 1.


Figure 1

Warning

- a) Install or remove the contra angle only after the motor has completely stopped running.
- b) Connect only to the adapter motor with a maximum speed of less than 40,000 rpm.

c) The connector of this product conforms to ISO 3964-2016 standard.

2.2 Mounting/dismantling of the needle (Figure 2)

- a) Insert the needle until it is in the correct position;
 - b) Press the button on the head and insert the needle into the chuck until the needle snaps into the mechanical recess, release the button;
 - c) When not pressing the button, gently push and pull the needle to confirm that the needle is secure.
- When disassembling, press the button first, then pull out the needle.

Figure 2

Notice

- a) After the needle is locked, gently pull on the needle to ensure that it is locked.
- b) If the needle is not securely mounted or not in place, it may fly out or not be removed.

3 Cleaning, disinfection and sterilization

Steps	Reprocess of dental low speed handpiece
Warnings	1. Before use, please carefully read the Instruction Manual of handpiece. 2. Before first use, please must clean and sterilize the Contra-angle Handpieces. 3. The use of ultrasound cleaning device and strong cleaning and disinfection fluids (alkaline pH>9 or acid pH <5) can reduce the life span of products. The manufacturer takes no responsibility in such cases. 4. The products may not be exposed to temperature above 138°C. This device shall not be exposed to high temperature above 138°C.

	5. After the device is damaged or reaches the Processing limit, please cleaned and sterilized the Contra-angle Handpieces before scrapping.
Processing limit	The products have been designed for a large number of sterilization cycles. The materials used in manufacture were selected accordingly. However with every renewed preparation for use, thermal and chemical stresses will result in ageing of the products. The allowed maximum times of sterilization for handpiece is 600 times.
Initial processing	1. Processing principles It is only possible to carry out effective sterilization after the completion of effective cleaning and disinfection. Please ensure that, as part of your responsibility for the sterility of products during use, only sufficiently validated equipment and product-specific procedures are used for cleaning/disinfection and sterilization, and that the validated parameters are adhered to during every cycle.

Please also observe the applicable legal requirements in your country as well as the hygiene regulations of the hospital or clinic, especially with regard to the additional requirements for the inactivation of prions.

2. Post-operative treatment

The post-operative treatment must be carried out immediately, no later than 30 minutes after the completion of the operation. The steps are as follows:

- 1) Remove the contra angle from the motor handpiece, and rinse away the dirt on the surface of contra angle with pure water (or distilled water/deionized water);
 - 2) Dry the contra angle with a clean, soft cloth and place it in a clean tray.
- Cautions:

	- The water used here must be pure water, distilled water or deionized water. Use a clean soft brush to carefully brush the head and back cover of the contra-angle until the dirt on surface is not visible. Then use soft cloth to dry the contra-angle and accessories and put it into a clean tray. The cleaning agent can be pure water, distilled water or deionized water.
Transportation	Safe storage and transportation to the reprocessing area to avoid any damage and contamination to the environment.
Cleaning	Preparation before cleaning Steps: Tools: tray, soft brush, clean and dry soft cloth 1. Remove the shanks/files and disconnect the contra-angle from the motor handpiece, put the contra-angle into the tray (Refer to article 2: Operation instructions)

	2. Use a clean soft brush to carefully brush the head and back cover of the contra-angle until the dirt on surface is not visible. Then use soft cloth to dry the contra-angle and accessories and put it into a clean tray. The cleaning agent can be pure water, distilled water or deionized water. The cleaning should be performed no later than 24 hours after the operation (neodisher MediClean forte).
	The cleaning can be divided into automated cleaning and manual cleaning. Automated cleaning is preferred if conditions permit. Automated cleaning Use a washer-disinfector meeting the requirements of the ISO 15883 series.

Carefully place the contra-angle in the disinfection basket. Fastening of the contra-angle if only permissible if they are freely moveable in the fixture. The contra-angle are not permitted to make contact with one another. Use a suitable rinsing adaptor, and connect the internal water lines to the rinsing connection of the washer-disinfector and start the program:

- 4min pre-washing with cold water(<40°C);
- Emptying
- 5min washing with a mild alkaline cleaner at 55°C
- Emptying
- 3min neutralizing with warm water(>40°C);
- Emptying

	• 5min intermediate rinsing with warm water(>40°C); • Emptying The automated cleaning processes have been validated by using 0.5% neodisher MediClean forte(Dr. Weigert). Note Acc. To EN ISO 17664 no manual reprocessing methods are require for these devices. If a manual reprocessing method has to be used, please validate it prior to use. The used cleaner is 0.5% neodisher MediClean forte.The intrinsic suitability of the contra-angle for effective cleaning and disinfection using the above automated cleaning and disinfection procedures was verified by a certified facility. (The used cleaner is 0.5% neodisher MediClean forte).
Disinfection	Not applicable.

Drying

If your cleaning and disinfection process does not have an automatic drying function, dry it after cleaning and disinfection.

1. Spread a clean white paper (white cloth) on the flat table, point the contra angle against the white paper (white cloth), and then dry the contra angle with filtered dry compressed air (maximum pressure 3 bar). Until no liquid is sprayed onto the white paper (white cloth), the contra angle drying is completed.

2. It can also be dried directly in a medical drying cabinet (or oven). The recommended drying temperature is 80°C~120°C and the time should be 15~40 minutes. Precautions:

- 1) The drying of product must be performed in a clean place.
- 2) The drying temperature should not exceed 138 °C;
- 3) The equipment used should be inspected and maintained regularly.


Inspection and maintenance	In this chapter, we only check the appearance of the handpiece. 1. Check the handpiece. If there is still visible stain on the handpiece after cleaning/disinfection, the entire cleaning/ disinfection process must be repeated. 2. Check the handpiece. If it is obviously damaged, smashed, detached, corroded or bent, it must be scrapped and not allowed to continue to be used. 3. Check the handpiece. If the accessories are found to be damaged, please replace it before use. And the new accessories for replacement must be cleaned, disinfected and dried. 4. If the service time (number of times) of the handpiece reaches the specified service life (number of times), please replace it in time.
-----------------------------------	---

Injection

Spray lubricating oil onto the Contra-angle Handpieces after cleaning and drying: The nozzle of cleaning lubricant is aligned with the air intake hole at the end of the contra angle to inject oil for 1-2 seconds.

Dental low speed handpiece

nozzle

Packaging	The disinfected and dried contra-angle and their accessories are assembled and quickly packaged in a medical sterilization bag (or special holder, sterile box). Precautions 1) The package used conforms to ISO 11607; 2) It can withstand high temperature of 137 °C and has sufficient steam permeability; 3) The packaging environment and related tools must be cleaned regularly to ensure cleanliness and prevent the introduction of contaminants; 4) Avoid contact with parts of different metals when packaging.
-----------	---

Sterilization	Use only the following steam sterilization procedures (fractional pre-vacuum procedure*) for sterilization, and other sterilization procedures are prohibited: 1. The steam sterilizer complies with EN13060 or is certified according to EN 285 to comply with EN ISO 17665; 2. The highest sterilization temperature is 138 °C; 3. Sterilization parameters: 5 min at 134°C, drying time: 20 min. Verification of the fundamental suitability of the products for effective steam sterilization was provided by a verified testing laboratory. Precautions: 1) Only products that have been effectively cleaned and disinfected are allowed to be sterilized;
---------------	---

- | | |
|--|---|
| | <p>2) Before using the sterilizer for sterilization, read the Instruction Manual provided by the equipment manufacturer and follow the instructions.</p> <p>3) Do not use hot air sterilization and radiation sterilization as this may result in damage to the product;</p> <p>4) Please use the recommended sterilization procedures for sterilization. It is not recommended to sterilize with other sterilization procedures such as ethylene oxide, formaldehyde and low temperature plasma sterilization. The manufacturer assumes no responsibility for the procedures that have not been recommended. If you use the sterilization procedures that have not been recommended, please adhere to related effective standards and verify the suitability and effectiveness.</p> <p>*Fractional pre-vacuum procedure = steam sterilization with repetitive pre-vacuum. The procedure used here is to perform steam sterilization through three pre-vacuums.</p> |
|--|---|

Precautions	1) Only products that have been effectively cleaned is allowed to be sterilized; 2) Before using the sterilizer for sterilization, read the Instruction Manual provided by the equipment manufacturer and follow the instructions. 3) Do not use hot air sterilization and radiation sterilization as this may result in damage to the product; 4) Please use the recommended sterilization procedures for sterilization. It is not recommended to sterilize with other sterilization procedures such as ethylene oxide, formaldehyde and low temperature plasma sterilization. The manufacturer assumes no responsibility for the procedures that have not been recommended. If you use the sterilization procedures that have not been recommended, please adhere to related effective standards and verify the suitability and effectiveness.
-------------	--

	*Fractional pre-vacuum procedure = steam sterilization with repetitive pre-vacuum. The procedure used here is to perform steam sterilization through three pre-vacuums.
Storage	Storage of sterilized instruments in a dry, clean and dust free environment at modest temperatures, refer to label and instructions for use.
Reprocessing validation study information	It is the duty of the user to ensure that the reprocessing processes including resources, materials and personnel are capable to reach the required results. State of the art and often national law requiring these processes and included resources to be validated and maintained properly.

4 Maintenance

When not in use for a long time, the following operations should be carried out once a week: oil lubrication, then manually check whether the needle can rotate normally; work once every 2 minutes at no load;

5 Storage and transport conditions

5.1 Storage and transport condition

5.1.1 Environment temperature: -20° C~+55° C

5.1.2 Relative humidity: 10%~93%

5.1.3 Atmospheric pressure: 70kPa~106kPa

5.2 This product should be placed or stored in a dry and clean environment which is away from acid, alkali and other harmful chemicals and gases.

5.3 Sterilized instruments should be kept in sterilization bags and kept intact and sealed. If stored for more than 7 days or if the sealed bags are damaged, the instruments should be reprocessed.

6. Environment protection

Patr	Toxic or harmful substances or elements					
	Pb	Hg	Cr6	Cr6+	PBB	PBDE
Main unit	o	o	o	o	o	o

1. The Contra-angle Handpieces meets EU RoHS environmental protection requirements.
2. When the handpiece is damaged or reaches the limit of use, the operator should clean and sterilize it in strict accordance with the requirements for reprocessing, and then hand it over to the local dealer or manufacturer;
3. Local distributors or manufacturers should scrap and recycle handpiece as metal products when they receive damaged or reach the limit of use.

7 After-sales service

After the sale of the product, if there is any quality problem and the device cannot work normally during warranty period, Woodpecker is responsible for the maintenance. Please refer to warranty card for the warranty period and other specific items.

This product is a precision equipment. If there is problem that needs to be repaired, returned to Woodpecker or handled by professionals is recommended.

8 Symble

Manufacturer

Type B applied part

Date of manufacture

Follow Instructions for Use

Used indoors only

Serial number

Fragile items, handle with care

Recyclable

IPX0 General equipment

Products comply with WEEE directive


Autoclavable up to the specified temperature

Avoid exposure to rain

Humidity limit for storage: 10% ~ 93%

Temperature limit for storage:
-20° C ~ +55° C

Medical Device

Keep away from sunlight

Atmospheric pressure for
operation and storage: 70kPa— 106kPa

CE marked product

EU Importer

Catalogue number

Authorised Representative in the
EUROPEAN COMMUNITY

9 Statement

9.1 If there are any issues not covered (such as: a notice to the user and/or patient that any series incident), the woodpecker will be announced on the website www.glwoodpecker.com

Guilin Woodpecker Medical Instrument Co.,Ltd.

Information Industrial Park, Guilin National High-Tech

Zone, Guilin, Guangxi, 541004 P. R. China

Sales Dept.: +86-773-5873196/2350599

After-sales Service Dept.: 0773-5827898

E-mail: woodpecker4@glwoodpecker.com

Website: <http://www.glwoodpecker.com>

MedNet EC-REP C IIb GmbH

Borkstrasse 10 · 48163 Muenster · Germany

ZMN-SM-1019 V1.0-20250210

PROCLINIC S.A.U.

C/ Palermo 9

50197 Zaragoza (España)

Rev. 02/25