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tray, guide, I-B-T 385, gingiva 2.0, cast NF,
match, model, model align, sprint, express

EN - Instructions for use

1. Product description

optiprint products are light-curing 3D printing resins for the production of dental molded parts (see second intended purpose) as a custom-made product for various (medical) indications, prescribed by dental professionals. They are suitable for use in devices of 385 nm and 405 nm wavelength unless otherwise stated on the label. The target patient group of the medical devices is defined as adults and adolescents.

2. Intended purpose

3D printing resin for the additive manufacturing of:

optiprint tray / optiprint tray 405	Individual impression trays		
optiprint guide 385 / optiprint guide HR / optiprint guide 405	Drilling templates		
optiprint tray viper green optiprint I-B-T 385	Transfer templates for bracket positioning		
optiprint gingiva 2.0	Gingival masks		
optiprint cast NF	Dental cast objects		
optiprint match	Scanable dental molded parts		
optiprint model / optiprint sprint / optiprint express	Dental models		
optiprint model align	Aligner Models		

The range also includes: optiprint clara, optiprint lumina and optiprint laviva.
Separate instructions for use are available for these products and must be observed.

3. Contraindications

The material should not be used for any purpose other than additive manufacturing of the intended purpose. Do not use the polymer material in the case of allergies to any of the ingredients (contains methacrylate monomers and oligomers).

4. Risk minimization and safety instructions

• Improper use and deviations from the described processing will result in impaired quality and biocompatibility as well as undesirable mechanical properties of the finished molded part.

• Biocompatibility is only guaranteed with proper application (post-polymerization without oxygen). Only process all molded parts when they are fully polymerized.

• The post-polymerization of the molded parts takes place in a suitable polymerization device (e.g., otofash G171), see 5.4.

• After completion of the construction process, the molded part should be cleaned with suitable cleaning solution (e.g., isopropanol 99 %) in an ultrasonic tank.

• The batch no. must be indicated in any procedure requiring the material to be identified.

• For additively manufactured medical devices made of printed resins, storing in water for 24 hours is recommended.

• Observe all recommended settings of the printer and the light-curing device.

• Read and observe the Safety Data Sheet (SDS) before use.

• For rupture safety, usual material wall thicknesses must be observed when creating the custom-made products.

• Proper personal protective equipment (nitrile gloves, safety glasses, protective clothing) must be worn when handling optiprint products as well as uncured printed molded parts.

• Avoid contact with skin and eyes before post-curing. The optiprint product can irritate eyes and skin.

• In rare cases, allergic reactions to components of optiprint products may occur. In case of accidental contact, follow first aid measures (rinse thoroughly with water and consult a doctor if necessary). See the SDS.

• To be used by qualified professionals only. Keep out of reach of children!

5. Processing steps

5.1 Designing

Please observe general design rules when designing molded parts. For optiprint tray, we recommend a minimum wall thickness of 2.5 mm. For hollow dental models, we recommend a minimum wall thickness of 3.5 mm.

5.2 Printing

Observe the instructions for use of the printer and the software. Ensure that work is carried out in clean conditions. Contamination on the 3D printer can cause defects on the molded part and damage the tray. The relevant material parameters for optiprint printing resins can be downloaded from the dentona homepage. Shake the optiprint bottle before use, and fill the material tray of the printer with sufficient material. Remove any bubbles with a clean instrument/spatula.

If you are using a printer without integrated heating, it is recommended to bring both the printer and the optiprint product to an operating temperature of 30 °C /86 °F. A cold start should be avoided.

5.3 Cleaning after printing:

Strip off any residue of the printing resin with a clean wooden spatula. Then remove the building platform from the printer and detach the molded parts from the building platform. Disconnect the support structures before cleaning.

Recommendation: In a non-heated ultrasonic bath: place cleaning fluid and the molded part in a second, sealable container (2 x 5 minutes – cast NF 2 minutes). Isopropanol (99 %) and optiprint clean can be used as cleaning fluids. When using optiprint clean, a 2-minute follow-up cleaning with isopropanol is always required. optiprint clean is not suitable for cleaning surfaces and equipment. Following this, dry the molded part with compressed air and check the results of the cleaning.

Complete the post-treatment of the molded parts quickly and adhere to the time specifications! In particular, avoid an unnecessarily long bath in the cleaning fluid and long waiting times between steps!

Alternatively, use the corresponding Rapid Shape Wash program with isopropanol (99 %).

You can determine cleaning has been sufficient if the surface is matt. Shiny areas require spot post-cleaning with cleaning fluid and a brush. Clean until no shiny areas are visible.

5.4. Finishing off

The final properties and also the final color depend on the post-curing process. To achieve the desired material properties and biocompatibility as well as the coloration, the fully cleaned and dried molded parts must be post-cured in oxygen-free conditions.

Recommendation

Otofash G171 flash-light device (NK Optik) with nitrogen flooding.

Ensure you use the Plexiglas tub with a UVB filter! It can be identified by the inscription: NK Optik 360N2.

Alternatively, use the corresponding Rapid Shape Cure program with a vacuum.

When post-curing in other light curing devices, make sure that the curing device has a product-specific, predefined setting or a comparably high energy input (200 W).

OTOFASH G171 LIGHT OUTPUT

The object is to be turned in between.

optiprint tray, guide, I-B-T 385, gingiva 2.0, match, model, model align, express	2 x 2000 flashes
optiprint tray viper green optiprint sprint	2 x 3000 flashes
optiprint cast NF	2 x 500 flashes
Alternative light output (200 W, 315-400 nm)	7 minutes (cast NF 2 minutes)

6. Polishing and cleaning

Dental molded parts produced with optiprint can be polished in the conventional way.

Final cleaning

A final cleaning is performed with detergent and water in a cold ultrasonic bath (5 minutes) until no polishing agent residues can be seen and felt on the molded part. Then rinse the molded part in clear water for 30 seconds until no more foaming from the detergent is visible.

7. Instructions for the laboratory/dentist for delivery of the molded part to the patient

Storing in water for 24 hours is recommended before inserting the molded parts. Disinfection is not necessary as standard.

8. Technical data

Please see the separate TDS (Technical Data Sheet).

9. Composition

Mixture of acrylate and methacrylate resins, photoinitiators (phosphine oxide), additives, dyes.

10. Storage

Protect the product from strong light and heat sources, recommended storage temperature: 5 °C to 30 °C. Close the bottle after each use. Do not use after the expiry date.

11. Disposal

Dispose of liquid printing resin in accordance with official regulations. Must not be disposed of with household waste. Do not discharge into drains.

12. Information

Report all serious incidents (death, serious deterioration of health, serious risk to public health) related to the medical device to the manufacturer and the competent authority of the Member State.

If defects in the material occur within the warranty period, the user is only entitled to replacement of the material if the defect is the responsibility of the manufacturer. dentona AG shall not be liable for any loss or damage caused by this material, whether direct or indirect, or special incidental or consequential damages, regardless of the legal basis. dentona AG is liable for direct material damage to the material caused by intent or gross negligence on the part of its legal representatives or executive employees, as well as for personal injury in accordance with the statutory provisions. Any liability for the material and consequential damages resulting from its use are excluded if the user has not followed the specified steps of the procedure.

13. Symbols

 Please follow the instructions for use (also available at www.dentona.de)

 Use by/expiry date

 Protect from sunlight

 Caution: hazardous to health

 Caution: corrosive

 Medical device

 Batch number

 Unique Formula Identifier

 Article number

 Caution: systemic health hazard

 CE marking

 Manufacturer

 Temperature restriction (5-30 °C)

 Caution: hazardous to the environment

 Unique Device Identification

tray, guide, I-B-T 385, gingiva 2.0, cast NF,
match, model, model align, sprint, express

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