vivapink

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Chargennummer

Indication:

Dental resin for the fabrication of all kinds of flexible dental therapeutic interim restorations with clamps. For repair, relining and fixation of dentures in denture bases made of dentona® vivapink, optimill memopink or optiprint interim.

The use of the interim restoration is limited to a maximum of 30 days.

Contraindication:

The product should not be used for any other purpose than the fabrication, repair or relining of (clamp) interim restorations and the fixation of acrylic dentures to these. dentona® vivapink is not suitable for the fixation of ceramic dentures. The product is not compatible with hard denture resins (thermosets). Unwanted biological reactions (such as allergies against material components) can occur in very rare cases. At known incompatibilities or known allergies to the material components works made from dentona® vivapink must not be used.

Hazard warnings:

The use of dentona® vivapink can produce harmful dusts that can irritate the eyes, skin and airways. Always ensure that the extraction system at your workplace is working properly (see also safety data sheet at www.dentona.de).

Processing:

Mixing ratio for dentona® vivapink: 1 part powder / 0.75 parts liquid

A complete restoration requires 10 g powder / 7.5 g liquid.

The mixing ratio and the temperature during processing (5-9 °C) must be strictly observed. Both components should be refrigerated for at least 15 minutes before mixing. Both components should be shaken for 60 seconds before mixing. Only use mixing cups and mixing spatulas that do not come into contact with any other resin.

Pack pressing technique:

1. Prepare the model in the normal way, measure and mark the tooth equator. As dentona® vivapink is a flexible material, the work can be arranged slightly above the tooth equator.
2. The work to be created is waxed-up in the articulator (for example with pink wax sheet) taking into account the chewing movement, and then bed the model into a cuvette. Now cool the required amount of powder in a refrigerator. Both components must be cooled before mixing.
3. Isolate all exposed plaster parts of the embedding with alginate isolation. As soon as the insulation is dry the counterpart can be produced from the modelled wax work. Once the counterpart gypsum has completely set, the entire flask should be warmed up in a warm water bath. This makes it easier to open the two halves of the cuvette.
4. Remove the wax from the modelled wax work, scald the model and the counterpart, water and isolate with alginate insulation.

IMPORTANT: Adhere exactly!

5. Remove the dentona® vivapink components from the refrigerator immediately prior to processing.
6. Shake both dentona® vivapink components for 60 seconds.
7. Mixing ratio: 1 part powder / 0.75 parts liquid.
8. Measure the dentona® vivapink material in separate dosing cups. Bring together in a rubber mixing cup, premix rapidly and bubble-free for 5-10 seconds and leave covered.
9. As soon as dentona® vivapink has a pasty consistency after approx. 5 minutes, the resin is placed in the cuvette on the row of teeth and the cuvette top is attached. Any visible white dots in the mixed paste are production-related and will no longer be visible after polymerization.
10. Now quickly screw the cuvette under the press and then polymerize for at least 45 minutes in a water bath at 2 bar and 50 °C (maximum 53 °C).
11. The cuvette must be completely cooled down before opening. This is ensured by placing the cuvette in a cold water bath for approximately 15 minutes.
12. Now carefully lift the work from the model and perform the usual elaborating and polishing.

Casting technique:

1. Prepare the model in the normal way, measure and mark the tooth equator. As dentona® vivapink is a flexible material, the work can be arranged slightly above the tooth equator.
2. The work to be created is waxed-up in the articulator (for example with pink wax sheet) taking into account the chewing movement. Now cool the required amount of powder in a refrigerator. Both components must be cooled before mixing.
3. Now prepare a silicone counterpart (for example dentona® 1:1 gum) and allow to cure. Once the silicone counterpart is fully cured it can be carefully lifted off. Cut out an excess channel with a sharp scalpel on the tuber or both sides on the retromolar triangle.
4. Remove the wax from the modelled wax work, scald the model and the counterpart, water and isolate the still warm plaster model with alginate insulation. The silicone counterpart is then best isolated with a thin layer of vaseline. Once the insulation has dried and the model has completely cooled down (simply place it to the resin in the refrigerator), place the counterpart back on the model and fix securely.

IMPORTANT: Adhere exactly!

5. Remove the dentona® vivapink components from the refrigerator immediately prior to processing.
6. Shake both dentona® vivapink components for 60 seconds.
7. Mixing ratio: 1 part powder / 0.75 parts liquid.
8. Measure the dentona® vivapink material in separate dosing cups. Bring together in a rubber mixing cup, premix rapidly and bubble-free for 5-10 seconds and then immediately fill into one side of the counterpart until the resin gets visible on the other side.
9. Now polymerise the model with the filled silicone counterpart as quickly as possible for 45 minutes in a pressure pot at 2 bar and a minimum of 50 °C (maximum 53 °C).
10. Then the work is removed from the pressure pot and cooled in cold water before the finished work is taken off. Now carefully lift the work from the model and perform the usual finishing and polishing.

Fixation of dentures:

1. The corresponding area of the denture base must be roughened and wetted with a small amount of dentona® vivapink liquid.
2. Undercuts should be inserted into the dentures.

IMPORTANT: Adhere exactly!

3. Remove the dentona® vivapink components from the refrigerator immediately prior to processing.
4. Shake both dentona® vivapink components for 60 seconds.
5. Mixing ratio: 1 part powder / 0.75 parts liquid
6. Measure the dentona® vivatemp material in separate dosing cups. Bring together in a rubber mixing cup, premix rapidly and bubble-free for 40-60 seconds. If necessary, leave to stand covered until the desired consistency is reached for further processing.
7. Polymerise for 45 minutes in a pressure pot at 2 bar and a minimum of 50 °C (maximum 53 °C).
8. Remove the work from the pressure pot and cool in cold water.
9. Elaborate and polish the work as usual.

Repair and relining:

1. The corresponding area of the denture base must be roughened and wetted with a small amount of dentona® vivapink liquid.
2. Soak the model in water and isolate the still warm plaster model with alginate insulation; it is best to isolate any silicone counters with a thin layer of vaseline. Once the insulation has dried and the model has completely cooled down, simply place it to the resin in the refrigerator.

IMPORTANT: Adhere exactly!

3. Remove the dentona® vivapink components from the refrigerator immediately prior to processing.
4. Shake both dentona® vivapink components for 60 seconds.
5. Mixing ratio: 1 part powder / 0.75 parts liquid
6. Measure the dentona® vivatemp material in separate dosing cups. Bring together in a rubber mixing cup, premix rapidly and bubble-free for 40-60 seconds. If necessary, leave to stand covered until the desired consistency is reached for further processing.
7. Polymerise for 45 minutes in a pressure pot at 2 bar and a minimum of 50 °C (maximum 53 °C).
8. Remove the work from the pressure pot and cool in cold water.
9. Elaborate and polish the work as usual.

Application tips:

1. Before removing from the model and before inserting in the patient's mouth the work should be made slightly flexible by exposing to lukewarm water.
2. For pre-polishing, use a leather buff with pumice at 3000 rpm. For the high-gloss polish, the cotton swab with a conventional high-gloss polish at 1500 rpm is suitable.
3. Severe undercut areas and interdental spaces should be sufficiently blocked out so that after finishing the restoration can be lifted from the model without damage.
4. The minimum height of the clasp arms after finishing must not be less than 2.0 mm!
5. The minimum thickness of the clasp arms after finishing must not be less than 1.0 mm!

Care/Cleaning of the restoration:

The cleaning of the restoration should not be done by excessive evaporation due to the high heat development. Clean the objects for approx. 5 minutes in water in an ultrasonic bath at a temperature of 35 °C until no milling or polishing residue can be seen or felt.

Cleaning solutions, alcoholic desinfectants and denture cleaners are not recommended. Wearing the restoration while drinking alcoholic beverages is not recommended. After meals, the restoration should be rinsed under cold running water.

Notes:

Store in clear water for at least 24 hours before delivery. Specify LOT No. for any operation that requires identification of the material.

Composition:

Mono and polymethacrylate, stabilizer, radical initiator, catalyst

Storage temperature before opening: 5-30 °C

Storage recommendation for the restoration: Avoid direct sunlight.

Disposal:

Dispose of dentona® vivapink in accordance with official regulations. Must not be disposed of with household waste. Do not discharge into drains.

Use:

Use only by dental professionals. Rx only.

Patient target group: adults and adolescents

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