Icon

Inhalt | Contents

DE – Gebrauchsinformation	4
EN – Instructions for use	14
FR – Mode d'emploi	23
IT – Istruzioni d'uso	34
ES – Instrucción de uso	44
PT – Instruções de uso	54
NL – Gebruiksaanwijzing	65
DA – Brugsanvisning	75
SV – Användarinstruktioner	84
PL – Instrukcja użycia	94
RU – Информация по применению	104
ZH – 使用说明	115
KO – 사용 정보	123

Produktbeschreibung

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Product description

Icon is an innovative product for the micro-invasive treatment of cariogenic lesions in proximal regions and on smooth surfaces and available in two variations:

Icon Caries Infiltrant – Proximal is a specially developed hard tissue preserving treatment for incipient proximal caries.

Icon Caries Infiltrant – Smooth Surface is a specially developed infiltration product for smooth surfaces and particularly well-suited to remove white spots and treat incipient caries for orthodontic patients after bracket removal.

Indications

- Micro-invasive treatment of non-cavitated cariogenic enamel lesions (lesion depth up to D1*).
- Esthetic improvement of mild to moderate non-cavitated fluorotic lesions (on smooth surfaces)

Lesion depth classification*

* Radiographic lesion depth classification according to bitewing images.

Contraindications

Do not use this material:

- for deeper seated lesions (D2 – D3) (see »Lesion depth classification*« on page 14) with cavitated enamel (enamel defects).
- in case of known allergies to any material component or existing contact allergies.

General safety information

- ▶ For dental use only!
- ▶ Keep out of reach of children!
- ▶ **Caution!** Icon-Etch can cause chemical burns. Avoid contact with soft oral tissue, eyes, and skin. In case of inadvertent contact rinse with copious amounts of water, and, if necessary, consult a physician.

- ▶ Limit the application to the surface to be treated. Avoid contact of the etching gel with the gingiva and adjacent teeth.
- ▶ Use only the syringe and application tip provided to apply the etching gel in the patient's mouth.
- ▶ Discard the application tip after use and cap the syringe.

Caution! Avoid contact of Icon-Infiltrant and Icon-Dry with eyes and skin.

- ▶ In case of eye contact, rinse with copious amounts of water, and, if necessary, consult a physician.
- ▶ In case of skin contact, wash the areas thoroughly with water and soap.

Caution! Improper use or incorrect application of Icon-Etch may cause pain in patients.

- ▶ Use Icon-Etch on enamel surfaces only. Do not apply onto exposed dentine or root cementum.

Please note

- Light-curing units should have an output of 450nm and should be checked regularly.
- The light intensity should be at least 800mW/cm².
- Place the light unit as close to the material as possible.

- The etched surface must not be touched or contaminated with saliva or moisture until the treatment resumes. If contamination occurs after drying, re-etch for approx. 10 s and repeat drying with Icon-Dry.
- For best treatment results a sufficiently dry working area is critical. Therefore, appropriate measures for a complete isolation must be taken. Do not use rubber dams made from thermoplastic elastomers such as Flexi-Dam (Coltène Whaledent/Hygenic).
- Discard syringes after use.
- On superficially remineralized, mostly older white spots, the pseudo-intact surface layer may be fairly thick. In this case, a 2 minute etching step to remove the surface layer may not be sufficient and etching and rinsing may have to be repeated. The teeth are then dried with Icon-Dry and infiltrated.

Proximal surface treatment

- The syringes included in the treatment unit contain sufficient material for treating two proximal lesions.
- Additional proximal surfaces can be treated consecutively during one visit.

Smooth surface treatment

- The syringes included in the treatment unit contain sufficient material for treating two to three smooth surface lesions.

- Multiple smooth surface lesions can be treated at the same time.
- For smooth surface applications a liquid dam (Liquid Dam) can be used as well.
- Treating orthodontic patients after bracket removal: Unless white spots are treated early, i. e. 1 to 2 months after bracket removal, it is recommended to repeat the etching step. A third etching process is recommended if a white spot is still visible after Icon-Dry is applied.
- In total, a lesion can be etched up to three times for 2 minutes each with Icon-Etch.
- For topical cariogenic white spots an area of 2 mm beyond the lesion site is etched. At the discretion of the dentist the entire smooth surface area ought to be etched and infiltrated in case of large area white spots as they occur after bracket removal.


Recommended use

Proximal surface treatment with Icon Caries Infiltrant – Proximal

Preparation

1. Before the start of treatment, clean the affected tooth and adjacent teeth. Remove any cleaning residue with water spray.
2. A sufficiently dry working area is critical. Therefore, appropriate measures like rubber dam or liquid dam must be taken. Please follow the manufacturer's specifications and instructions.


Separation

3. Introduce one of the enclosed dental wedges into the interdendum.
 - ▶ For better access to the proximal region the wedge handle can be bent or detached with a twisting motion.

In order to obtain sufficient separation of the teeth the wedges must be seated at an adequate depth in the interdentalium.

- ▶ A separation of approx. 50 µm is recommended to ensure a successful treatment. This can be done gradually: Insert the wedge to the point of resistance. Wait 3 to 5 s and gradually move the wedge deeper until sufficient separation is achieved.
- ▶ Leave the wedge in the proximal space during the entire treatment procedure.
- ▶ Alternatively, especially at very narrow proximal spaces, it is possible to separate the teeth with an orthodontic rubber ring.


Uncovering the lesion body

For best treatment results the hyper mineralized surface layer must be removed.


4. Screw the Proximal-Tip onto the Icon-Etch syringe, and introduce the application tip into the interden-

tium. Be sure to align the green side of the Proximal-Tip with the area to be treated.

- ▶ **The material can only be dispensed on the green side of the application tip.**


5. Apply an ample amount of Icon-Etch onto the lesion site (1½ to 2 turns of the shaft). Let Icon-Etch set for 2 minutes. The etching gel will be activated by slightly moving the applicator. Remove excess material.

6. Remove the application tip from the interdendum. Suction off Icon-Etch and rinse with water for a minimum of 30 s. Dry carefully with oil-free and water-free air.

Drying


For best treatment results it is necessary to dry the lesion in-depth.


7. Screw the application tip onto the Icon-Dry syringe, apply an ample amount of material onto the lesion and allow to set for 30 s. Dry thoroughly with oil-free and water-free air.

Infiltration

The lesion can be penetrated in-depth by the Infiltrant after complete drying.

- ▶ **Do not apply Icon-Infiltrant under direct operating light as this may cause the material to set prematurely.**
8. Screw a new Proximal-Tip onto the Icon-Infiltrant syringe, and introduce the application tip into the interdendum. Be sure to align the green side of the Proximal-Tip with the area to be treated.
- ▶ **The material can only be dispensed on the green side of the Proximal-Tip.**
9. Apply an ample amount of Icon-Infiltrant onto the lesion site (1½ to 2 turns of the shaft).
10. Let Icon-Infiltrant set for 3 minutes. The Infiltrant will be activated by slightly moving the applicator. Add material if necessary.
11. Remove the application tip from the interdendum. Remove excess material with dental floss.

12. Light-cure Icon-Infiltrant from all sides for at least 40 s (total).
13. Screw a new Proximal-Tip onto the Icon-Infiltrant syringe, repeat the application, and allow to set for 1 minute. Remove the application tip from the interdentium. Remove excess material with dental floss. Then light-cure from all sides for at least 40 s (total).
14. If a second lesion needs to be treated repeat steps 1 through 13 for the affected site.
15. Remove the wedge and rubber dam. Use polishing strips for the surface finish.

Smooth surface treatment with Icon Caries Infiltrant – Smooth Surface

Preparation


1. Before the start of treatment, clean the affected tooth and adjacent teeth. Remove any cleaning residue with water spray.

2. For best treatment results a sufficiently dry working area is critical. Therefore, appropriate measures like rubber dam or liquid dam must be taken. Please follow the manufacturer's specifications and instructions. Ensure that the treated lesion is completely accessible.

Uncovering the lesion body


For best treatment results the hyper mineralized surface layer must be removed.

3. Screw the Smooth Surface-Tip onto the Icon-Etch syringe.

4. Apply an ample amount of Icon-Etch onto the lesion site by turning the syringe shaft carefully, allow to set for 2 minutes and activate it by moving. Remove excess material with a cotton wad.

Note: Unless white spots are treated early, i. e. 1 to 2 months after bracket removal, it is recommended to repeat the etching step up to three times.


5. Suction off Icon-Etch and rinse with water for at least 30 s. Dry with oil-free and water-free air.

Visual Check

For best treatment results it is necessary to dry the lesion in-depth.

- ▶ **In this step a preview of the final result is shown. When wetted with Icon-Dry, the whitish-opaque coloration on the etched enamel should diminish. If this is not the case repeat the etching step once or twice for 2 minutes each, and rinse and dry the teeth again (steps 3-5).**

6. Screw the application tip onto the Icon-Dry syringe, apply an ample amount of material onto the lesion, and allow to set for 30 s.

Drying


For best treatment results it is necessary to dry the lesion in-depth.

7. After the visual check dry the lesion thoroughly with oil free and water-free air.


Infiltration

The lesion can be penetrated in-depth by the Infiltrant after complete drying.

8. Screw a new Smooth Surface-Tip onto the Icon-Infiltrant syringe.
- ▶ **Do not apply Icon-Infiltrant under direct operating light as this may cause the material to set prematurely.**

9. Apply an ample amount of Icon-Infiltrant onto the etched surface by turning the shaft.
10. Allow Icon-Infiltrant to set for 3 minutes. The Infiltrant will be activated by slightly moving the applicator. Add material if necessary.
 - ▶ **In case of deeper and larger defects the esthetic result can be improved by extending the exposure time up to 6 minutes.**
11. Remove excess material with a cotton wad and dental floss.

12. Light-cure Icon-Infiltrant for 40 s.
13. Screw a new Smooth Surface-Tip onto the Icon-Infiltrant syringe, repeat the application, and allow to set for 1 minute. Remove excess material with a cotton wad and dental floss, and light-cure for a minimum of 40 s.
14. Remove the rubber dam. Use polishing cups (or similar) for the surface finish.

Side effects

- In some cases, contact allergies with similar composite products have been reported.
- Contact of Icon-Etch with the oral mucosa causes a white surface coloration. This will subside after a few days.

Interaction

Polymerization inhibiting substances such as products containing eugenol (e. g. temporary cements) must not be used in combination with this material.

Composition

- Icon-Etch: Hydrochloric acid, pyrogenic silicic acid, surface-active substances

- Icon-Dry: 99% ethanol
- Icon-Infiltrant: Methacrylate-based resin matrix, initiators, additives


Storage

- ▶ Store at 2 to 25 °C / 36 to 77 °F!
- ▶ Protect from direct sunlight!
- ▶ Do not use after the expiration date!

Patient card + sticker

This material is not radiopaque. To document the infiltrated surfaces and lesion depths the enclosed patient card can be filled out at the time of treatment and follow-up visits and handed to the patient. In addition, the enclosed sticker supports the documentation in the patient's record.

1. Front: Enter the name of the patient and office stamp.
2. Tooth diagram: Mark the treated area.
3. Record the treated tooth [1], the treated surface [2] (d = distal, m = mesial, v = vestibular (buccal, cervical, labial), l = lingual or palatal), lesion depth*, and the treatment date [3].

4. Follow-up visit [4]: Record the lesion depth* and date of the follow-up appointment. Annual check-ups are recommended.

see »Lesion depth classification« on page 14.

Icon Caries Infiltrant – Proximal	
Starter Kit @ 2 treatment units	REF 220341
Package @ 7 treatment units	REF 220237

Icon Caries Infiltrant – Smooth Surface	
Starter Kit @ 2 treatment units	REF 220343
Package @ 7 treatment units	REF 220238

One treatment unit includes:

Icon Caries Infiltrant – Proximal

- 1 Icon-Etch syringe @ 0.3 ml fluid
- 1 Icon-Dry syringe @ 0.45 ml fluid
- 1 Icon-Infiltrant syringe @ 0.45 ml fluid
- Accessories

Icon Caries Infiltrant – Smooth Surface

- 1 Icon-Etch syringe @ 0.45 ml fluid
- 1 Icon-Dry syringe @ 0.45 ml fluid
- 1 Icon-Infiltrant syringe @ 0.45 ml fluid
- Accessories

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DMG Chemisch-Pharmazeutische Fabrik GmbH