HI-SEAL

Definition

Radiopaque, light-cure Cavity Liner

Safety notes

Warnings

HI-SEAL contains polymerizable methacrylate resins which may irritate skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible persons.

Avoid eye contact to prevent irritation and possible corneal damage. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.

Avoid skin contact to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and wash thoroughly with water and soap. In case of skin sensitization or rash, discontinue use and seek medical attention.

Avoid contact with oral soft tissues/mucosa to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If sensitization of mucosa persists, seek medical attention.

Precautions

Wear suitable protective eyewear, mask, clothing and gloves.

Protective eyewear is recommended for patients.

Contraindications

HI-SEAL is contraindicated in patients with a history of allergic reaction to methacrylate resins or any other of the components. Do not use HI-SEAL in cases of pulpar exposition, pulpitis.

Instruction for use

- Prepare the cavity.
- Rinse the cavity and slightly dry.
- Isolate treatment area from salivary contamination (cotton rolls, rubber dam etc...)
- Apply the material directly over the dentin through a disposable applicator. When increased bond strength to dentin is desirable, apply first your preferred dental adhesive to the dentin
- The first coating should not be more than 0.3 - 0.5 mm thick. Remove any material inadvertently applied over enamel margins.
- Cure for 20 seconds, keeping the tip of the LED light as close to the cured surfaces as possible.
- Apply a second layer, if indicated. The second layer can be up to 1.5mm thick. Cure for 30 seconds.
- Restore the tooth with your usual material. Hi-Seal is compatible with all commercially available dental restorative materials.

Note

If using HI-SEAL under water based cements such as glass ionomer, remove any uncured resin from the Hi-Seal surface by wiping it with a pledget of cotton moistened with alcohol.

Storage & shelf-life

Do not expose the product to sun-light or warm sources. Never store HI-SEAL close to products containing eugenol. Keep the products in a cool place, at a temperature between 12°C and 25°C. When stored under such conditions, HI-SEAL has a shelf-life of 24 months.

Presentations & references

HS-U: 4x1.5g syringes + 8 needle tips

Important notice

- Do not throw away the box and the instruction for use until the product is completely finished.
- High humidity could damage the packaging on which important instructions and legal information are printed.
- Storage at a too high temperature could cause to age prematurely and inhibit its polymerization.
- Storage at a too low temperature leads to the decomposition of the product which will render it completely unusable.

Disclaimer

These materials have been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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