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INSTRUCTIONS FOR USE

LIGHT CURE RESIN MODIFIED GLASSIONOMER LUTING CEMENT

Proclinc Expert Resin modified Glass Ionomer Cement is a light cured resin resin modified glassionomer luting cement. Beside its excellent fluoride release and the excellent biocompatibility as glassionomer cement Proclinc

Expert Resin modified Glass Ionomer Cement has also good chemical bonding to dentine and enamel and to the restoration and a tight seal at the dental margins. Because of its radiopacity it ensures easy postoperative diagnosis.

Indications

- Cementing of crowns, bridges, inlays and onlays (all types: metal, resin, ceramic fused to metal and ceramic)
- Liner under composite fillings

1. Tooth Preparation

Prepare tooth in usual manner.

Apply Proclinc Expert Resin modified Glass Ionomer Cement Primer with a suitable brush onto the moist enamel and dentine surfaces for 30 seconds with agitation. The material should build a homogeneous layer. Air thin gently for 10 seconds to remove the volatile components and to disperse the Proclinc Expert Primer. Then light cure with a suitable dental light unit for 10 seconds.

Pulp capping with Proclinc Expert Resin modified Glass Ionomer Cement is contraindicated. To deep areas of possible pulpal exposure apply a small amount of a calcium hydroxide liner.

2. Dispensing and Mixing

The powder/liquid ratio to achieve a suitable consistency is

1.7 / 1.0. This can be obtained by mixing one level scoop of powder and 3 drops of liquid.

For accurate dispensing of Proclinc Expert Resin modified Glass Ionomer Cement powder shake the bottle to loosen the powder. Overfill the spoon with the powder, level the powder for using the collar at the top of the bottle and carry it onto the mixing pad. Avoid compressing powder into the spoon with the inside wall of the bottle.

For dispensing of Proclinc Expert Resin modified Glass Ionomer Cement liquid turn the bottle vertically with the tip about 5 cm above the mixing pad. Steady your hand and squeeze the liquid gently to dispense one drop. If any bubbles are present, lightly tap the bottle with the fingers holding it. Discout drops that contain bubbles and are obviously not full-sized.

Use a small spatula to rapidly mix all the cement powder into the liquid. The mixed cement should be thixotropic and have a smooth consistency and glossy appearance. Total mixing time is 30 sec.

After use, tightly close both liquid and powder bottles to prevent exposure to moisture.

3. Cementing

Avoid water and saliva contamination during application and setting of the cement. Rubber dam is the best way of achieving tooth isolation.

Prepare the restoration according to manufacturer instructions.

Mix the required amount of Proclinc Expert Resin modified Glass Ionomer Cement and apply a thin layer (1 mm or less) to the prepared bonding surface of the restoration and seat immediately (working time 2:30 min from start of mixing at 23°C). Net setting time without any light is approx. 4:30 min.

Note: Higher temperatures will shorten the working time, lower temperatures will prolong the working time.

An overextended working time will cause the loss of adhesion to enamel and dentine.

Remove excess cement at the first setting stage.

Maintain isolation until the set of the cement is verified (approx. 4:30 min). Setting can be optimized with a 20 sec light cure.

4. Conclusive Notes

This product is to be applied only by a dental professional in the manner as described in this instruction.

Do not use Proclinc Expert Resin modified Glass Ionomer Cement with patients who show an allergy to the material. In case of allergic reactions immediately stop the application, and advise the patient to consult a physician. An operator, who has a history of allergy to glass-ionomer cements should not handle Proclinc Expert Resin modified Glass Ionomer Cement.

Do not allow the cement mixture to contact the oral tissues or skin. In case of contact, remove the material with absorbent cotton soaked in alcohol and rinse with water.

Avoid eye contact of the cement mixture. In case of contact, immediately flush with water and seek medical treatment.

Do not mix the powder or liquid of Proclinc Expert Resin modified Glass Ionomer Cement with any other glass-ionomer product.

Direct contact with eugenol-based products is to be avoided. Eugenol inhibits the setting of LC glass-ionomer cements.

5. Storage

Store Proclinc Expert Resin modified Glass Ionomer Cement in a cool and dark place at

4:25 °C (39:78 °F). Temperature should not exceed 25 °C (78 °F). Do not use after expiry date.

Warranty

PL Superior Dental Materials GmbH warrants this product will be free from defects in material and manufacture. PL Superior Dental Materials GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and PL Superior Dental Materials GmbH's sole obligation shall be repair or replacement of the PL Superior product.

Limitation of Liability

Except where prohibited by law, PL Superior Dental Materials GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

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