Fiche technique des HERO 642®

Technical datasheet HERO 642®

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Non-sterile: *Available in packs of 6 instruments.*
Available in assorted packs containing 3 instruments of each taper for no. 20 (L 21 mm), no. 025 (L 25 mm) and no. 030 (L 25 mm), i.e. in packs of 9 instruments.
Available in boxes of 10 packs.

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Geometric specifications:

The HERO 642® instruments have a .06, .04 and .02 taper.
They have a circular section with three cutting edges.
The HERO 642® instruments are available in a "Classics" version (handle: Ø 2.35) for use on a reducing endodontic contra-angle.
The instruments conform with the applicable paragraphs of the EN ISO 3630-1 standard.
The "Classics" handle conforms with the EN ISO 1797-1 standard.

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Specifications for use:

The HERO 642® instruments are used on contra-angles in continuous rotation.

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Maximum recommended number of uses under normal conditions:

Non-sterile instruments: 10 uses for all numbers.
If the file is not visually damaged.

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**Pre-disinfection, cleaning and sterilization:**

For hygienic and sanitary safety reasons, all root canal instruments marketed by MICRO-MEGA® must be pre-disinfected, cleaned and sterilized before each use on a patient (excepting sterile instruments).

This procedure is recommended for all instruments :

- before the first use for non sterile instruments*
- between each use on patient for re-usable instruments*

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Pre-disinfection: *Before the immersion of the instruments in the disinfectant bath sort them by type and material category in order to avoid electrolysis and corrosion. Comply with the instructions and concentrations recommended by the pre-disinfectant manufacturer (an excessive concentration can cause chemical damage to the materials). Respect the immersion time (do not leave the instruments in the disinfectant bath overnight). Only use disinfectant products recommended for the disinfection of dental root canal instruments (do not use products containing phenol, hydrogen peroxide, soda or sodium hypochlorite or extremely aggressive special disinfectants for burs).*

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Note: MICRO-MEGA® has tested the most common disinfectant products on the market. In view of the large number of available products, it is recommended to pre-test a product on a sample or to contact the manufacturer of the used product.

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Cleaning: *Brush the instruments for dentinal debris elimination. Do not clean the instruments loosely in an ultrasound device. Do not use soda-based solutions, solutions containing mercury salt or alkaline and acid solutions. The use of a washer/disinfector should be tested beforehand in order to check that the instruments are not deteriorated.*

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Rinsing: *Rinse the instruments abundantly using demineralized or distilled water to avoid corrosion and alterations during sterilization. Dry with dry air. Discard deformed instruments (broken, twisted, despiralized instruments or instruments showing marks conducive to breakage).*

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Wrapping : *Put the instruments in a sterilization bag compliant with EN ISO 11607 standard to be used in a steam sterilizer up to 134°C and 2.2 bar during 18 minutes. Follow the instruction of the sterilization bag's manufacturer.*

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Sterilization : *Sterilize the instruments wrapped in a sterilization bag. Use an autoclave (saturated steam sterilizer): at 134°C and 2.2 bar during 18 minutes.*

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Storage : *Check the integrity of the sterilization bag before using the instruments.*

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Normal function of the HERO 642® instruments:

The instruments are used for the shaping of the radicular root canal.

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Operating mode:


Start the treatment with the exploration of the root canal or an initial penetration into the root canal in order to identify the root canal trajectory and determine the working length. In case of difficult root canals, the passage of the HERO 642® with a .06 taper helps to improve the access in the coronal half.

Identify the level of root canal difficulty in order to choose the adapted HERO 642® sequence:

Simple root canals: straight with curves < 5°

Root canals of medium difficulty: curves > 10° and < 25°

Difficult root canals: curves > 25°

"Crown-down" method:

Conicité 6%

Simple root canals / blue sequence:

1) .06 taper

Preparation of the coronal access, identification of the root canal orifices and initial penetration applying conventional methods.

Adjust the (black) rubber stop on the HERO 642® instrument no. 030 with a .06 taper so that it corresponds to 1/2 or 2/3 of the working length.

Penetrate with the rotating instrument into the root canal, irrigate regularly and abundantly.

Perform a short up and down movement of low amplitude with normal drilling pressure.

Penetrate with subsequent pressure until 1/2 or 2/3 of the working length are reached

Conicité 4%

2) .04 taper

Adjust the (grey) rubber stop on the HERO 642® instrument no. 030 with a .04 taper so that it corresponds to the working length minus 2 mm.

Proceed in the same way as previously and combine the movement with a circumferential filing movement (surfacing).

If the instrument does not advance further down, perform a recapitulation with a control hand file in order to verify the attained penetration depth

Conicité 2%

3) .02 taper

Adjust the (white) rubber stop on the HERO 642® instrument no. 030 with a .02 taper so that it corresponds to the working length. Perform the same quick pumping movement until the working length indicated by the stop is reached without excessive pressure.

Combine with a circumferential filing movement.

Irrigate abundantly and dry prior to obturation.

A. Canal facile (droit - courbe < 5°)
SEQUENCE BLEUE

Diamètre	6%	4%	2%
No. 30	1/2 - 2/3 LT	LT moins 2 mm	LT

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Root canals of medium difficulty / red sequence:

The "crown-down" method remains the same, however, it is recommended to start with the instruments of diameter 25 while respecting the steps of the red sequence. Use of 5 instruments.

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Difficult root canals / yellow sequence:

The "crown-down" method remains the same, however, it is recommended to start with the instruments of diameter 20 while respecting the steps of the yellow sequence. Use of 6 instruments.

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Enlargement of the apical zone

It is possible to continue with .02 tapered instruments with diameter 035, 040 and 045 in the apical third if the root canal presents a morphology which so requires (necrosed root canal or open apex).

Diamètre	
Nr. 35	LT
Nr. 40	LT
Nr. 45	LT

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Disposal:

After use, instruments must be placed in a secure container, used to collect cutting or sticking instruments (like needles or disposable bistoury) as per good dentistry practices.

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Dispositif médical de classe IIA (Directive 93/42)

Medical device class IIA (93/42 Directive)

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VALIDATION :

FUNCTION : Regulatory Affairs Manager Micro-Mega

NAME : Philippe CONCHE

DATE (YYYY-MM-DD) : 2019-05-06

SIGNATURE :

FUNCTION : Project Manager Micro-Mega

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DATE (YYYY-MM-DD): 2019-05-06

SIGNATURE :