D_ Devices

Advanced systems for dental devices

D_BIOLOGICAL INDICATORS FOR STEAM STERILIZATION PROCESSES

EN_ User manual

ES_ Manual de instrucciones

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
figura 1-1

PRODUCT DESCRIPTION:

Self-Contained Biological Indicators for monitoring steam sterilization consist of:

A plastic pipe and cap. The cap with a hole is placed on top of the pipe to allow the sterilizing agent to penetrate into the indicator. Between pipe and cap there is a filter paper which is used to avoid pollution. A crushable media ampule which contains modified culture medium with a pH indicator. The modified culture medium will change from the initial Purple color to Yellow in the presence of bacterial growth.

An inoculated carrier of *Geobacillus stearothermophilus* ATCC 7953 with a population level of $\geq 10^5$ is located at the bottom closed end of the pipe.

The LBS144 / LBS288 Biological Indicator, when used with an incubator (incubation temperature 54-60°C), provides a test result after 24 hours(LBS144) or 48hours(LBS288) with an accuracy of $\geq 97\%$.

The label on the outside of the plastic pipe contains a Type 1 chemical indicator strip.

INDICATIONS FOR USE:

The LBS144 / LBS288 is utilized to monitor steam sterilization processes efficacy at 121°C to 135°C.

PERFORMANCE CHARACTERISTICS

Population	$\geq 1.0 \times 10^5$ per carrier
RIT	$\geq 97\%$
Purity	No evidence of contamination present in sufficient numbers to adversely affect the finished product.
Resistance	D value at 121 °C ± 1 °C ≥ 1.5 mins. Z value ≥ 6 °C.
Dimensions	10.9 mm(cap diameter) , 8.8 mm(pipe diameter), 46 mm(height)

STORAGE AND SHELF LIFE

- 10°C to 30°C, 30% to 80% Relative Humidity.
- Keep away from sunlight, Protect from sterilizing agent and radioactive sources, Do not dry.
- Shelf life: 36 months from the date of manufacture.

INDICATIONS FOR USE:

- Remove an appropriate number of LBS BIs from the box, reseal the box, and check each BI for:
 - Indication of a damaged ampule including low medium fill volume, wet or dried medium inside pipe, spore strip appearing wet or discolored.
 - Missing or damaged components including cap, cap filter, spore strip, medium ampule and plastic pipe. Dispose of any damaged or questionable units.
- Identify the BIs by labeling with process information, place the BI in a suitable test pack which is representative of the load. Place this test pack in the most challenging area of the sterilizer and run the sterilization program as usual.

- After sterilization is complete, remove from the sterilizer and allow the pack and biological indicator to cool for a sufficient time, at least 5 minutes.
- Retrieve the BIs from the test load.
- After exposure to steam, the chemical indicator strip changes from pink to black. The purpose of the chemical indicator strip is to distinguish whether the BI has been used; black does not indicate successful sterilization.
- Recheck each BI for:
 - Indication of a damaged ampule including low medium fill volume, wet or dried medium inside pipe, spore strip appearing wet or discolored. Dispose of any questionable units. Results obtained from damaged units should be considered suspect.
- Place the BI in an upright position, compress the plastic pipe in the crushing device to break the glass ampule. Confirm that the spore strip at the bottom of the BI is fully saturated with the culture medium; Be careful not to let the culture medium come into contact with the filter paper on the cap and incubate the indicator in the automatic reader. Use an unsterilized BI at least once per day or at every new batch of indicators and process it as above as a positive control tube.
- Incubate for 24 hours(LBS144) or 48hours(LBS288) to observe the biological culture results. If the color of the culture medium changes from purple to yellow, it indicates a positive result. If the color of the culture medium remains unchanged (purple), the BI is negative.

COMPLIANCE

- ISO 11138-1:2017 Sterilization of health care products-Biological indicators-Part 1: General requirements.
- ISO 11138-3:2017 Sterilization of healthcare products-Biological indicators-Part 3 Biological indicators for moist heat sterilization processes.
- ISO 11138-8:2021 Sterilization of healthcare products-Biological indicators-Part 8 Method for validation of a reduced incubation time for a biological indicator.
- ISO 11140-1:2014 Sterilization of health care products-Chemical indicators-Part 1: General requirements.

NOTE

- The Positive Control BI should show positive. If the Positive Control as does not show signs of growth, consider the test invalid.
- Short excursions outside the range of temperature recommended will not impact the performance of the LBS BIs.
- Do not use after the expiration date. Do not refrigerate.
- Autoclave for not less than 30 minutes at 121°C or per other validated disposal cycle prior to discard.
- If the incubation time exceeds 48 hours, ensure that the culture solution does not evaporate by sealing the holes on the cap

LISTER BIOMEDICAL CO.,LTD. Edificio 2, No.72 Shengl
Road, Yixiu, Anqing, Anhui, China. Tel. +86 556 5689070
Correo electrónico: info@listerbiomed.com www.listerbiomed.com

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LISTER BIOMEDICAL CO., LTD. Edificio 2, No.72 Shengl