# Soft relining material









---

## AUTOMIX SYSTEM

Instructions for use

---

- DE** Gebrauchsanweisung.
- EN** Instructions for use.
- ES** Instrucciones de uso.
- FR** Mode d'emploi.
- IT** Istruzioni per l'uso.
- PL** Instrukcja postępowania.
- PT** Instruções de uso.

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6 Min. 30 Sek.

Abbindezeit

## LONG-TERM SOFT RELINING MATERIAL FOR CHAIRSIDE USE, A-SILICONE BASED, COLD CURING

### 1. PREPARATION OF THE DENTURE

- Removal of existing relinings (Fig. 1): The old, existing relining has to be removed completely. Grind off the denture base approx. 1–2 mm, depending on the situation. Incorporate at the vestibular, lingual resp. dorsal borders a shoulder of approx. 2–3 mm in height with a depth of minimum 1 mm. Round off the borders and clean the denture thoroughly.
- Degreasing (Fig. 2): Degrease by means of pure medical alcohol and let dry.
- Applying the primer for adhesion (Fig. 3): Apply the primer uniformly and completely onto the grinded denture areas to which the material has to adhere. Avoid contact with eyes! Wear suitable eye protection! Let the primer dry for approx. 1 minute. After use close the bottle immediately. Saliva, humidity or monomer liquids prevent the adhesion of Proclinic Expert Soft Relining Material to the denture base.

### 2. USE OF DISPENSER/MIXING JET

Mount cartridge onto dispenser (is sold separately, reorder 0156) according to special instructions. Remove and dispose of closure cap, do not re-use! For exact flow control extrude slightly material until uniform flow from both orifices is achieved. Install mixing tip onto cartridge and extrude needed material. After each use leave the used mixing tip in place to seal the cartridge.

### 3. DIRECT CHAIRSIDE RELINING METHOD

- Coating the denture (Fig. 4): Prepare denture as described under point 1. Apply Proclinic Expert Soft Relining Material cartridge material uniformly by means of dispenser (see point 2) onto the vestibular, lingual resp. dorsal borders. Then coat completely the remaining denture with a layer of minimum 2 mm. All denture areas in contact with the gingiva have to be coated.
- Functional movements (Fig. 5): Seat the coated denture into the patient's mouth and let him close in normal occlusion. Setting is finished after approx. 5 minutes. Then remove the relined

denture and trim and polish as described under point 5.

### 4. INDIRECT METHOD IN THE LABORATORY

Invest denture to be relined into a flask with counter. Remove denture from the flask and prepare as described under point 1. Coat the model with an usual alginate based insulation. Extrude the required amount of Proclinic Expert Soft Relining Material from the cartridge by means of dispenser (see point 2) and apply it uniformly onto the denture area to be relined (the layer of Proclinic Expert Soft Relining Material should be approx. 1 mm more than the base material has been grinded off). Close flask with counter. Press with light pressure for approx. 30 minutes or place flask in a pressure pot at 40–45 °C/104–113 °F for approx. 10 minutes in the water bath. Then deflask the relined denture and trim and polish as described under point 5.

### 5. FINAL TRIMMING AND VARNISHING

- Smoothing of uneven spots: Remove excess material with a sharp scalpel or fine scissors (Fig. 6). Then smooth rough areas and borders by means of a grinding instrument with approx. 8000 r.p.m. under low pressure (Fig. 7).
- Polish (Fig. 8): Proclinic Expert Soft Relining Material polish is applied only for smoothing of the trimmed functional borders (mixing ratio: polish and catalyst 1:1). Avoid contact with eyes and skin! Wear suitable eye protection and gloves! Open the aluminium cap with a spatula, squeeze both liquids onto a mixing pad in a mixing ratio of 1:1, spatulate to a homogeneous mixture and brush on the mixed liquid. Let dry at room temperature for approx. 5 minutes.

### 6. CLEANING HINTS FOR THE PATIENT

Brush the relined denture several times a day under running water with a tooth brush and neutral soap. Cleaning bath: up to a maximum of 10 minutes.

## IMPORTANT WORKING HINTS

- Absolutely firm bonding is obtained only by applying the primer onto roughened surfaces.
- Primer must cover the entire area completely! For better application check the primer contains a colour indicator.
- Proclinic Expert Soft Relining Material coating must cover the surface of denture in contact with the gingiva completely and uniformly. Otherwise adhesion problems may occur.
- Setting times must be strictly observed.
- For trimming with grinding sleeve apply pressure with approx. 8000 r.p.m. – use grinding sleeve for silicone materials only.
- Please follow the instructions of the safety data sheets.
- Primer: avoid contact with eyes, mucous membranes, and skin. Do not inhale vapors.
- Polish: avoid contact of catalyst with eyes and skin, this can lead to irritation. In case of accidental contact, rinse eyes immediately with running water for at least 15 minutes, consult an ophthalmologist.

## INDICATIONS FOR USE

- Suited for all synthetics based on polymethylmethacrylate.
- Fast and easy fabrication of permanently soft relinings; chairside or indirect use in the dental laboratory.
- As elastic final margin of maxillary denture and to steady the position of the upper and lower dentures by increased adhesion.
- Post-damming of the A-line in case of adhesion problems and with flabby ridges.
- To steady the position of the upper and lower dentures by increased adhesion.
- Elastic cover of sharpridged alveolar processes.

- For covering implant regions and as overdentures.

## TECHNICAL DATA

DIN EN ISO 10139-2 Type A

- Mixed volume: 50 ml cartridges
- Colour code:  
Base: gingiva-coloured  
Catalyst: white
- Mixing ratio: 1:1
- Mixing time: does not apply (Automix-system)
- Working time: approx. 1½ min.\*
- Setting time: approx. 6½ min.\*
- Time in mouth: approx. 5 min.
- Linear dimensional change: < 0,3 %
- Recovery from deformation: > 99 %


- Application: At 23 °C/73 °F, 50 ± 5 % rel. humidity.
- Storage:  
\* from beginning of mixing at 23 °C / 73 °F, 50 ± 5 % rel. humidity. Increased temperatures accelerate, decreased temperatures retard the indicated times.

## Further information

Silicone based materials are proven a million times. On condition of a proper application, undesired effects are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to get tested the allergenic potential before the application of the material.

For use by trained specialist staff.

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**DETAX GmbH & Co. KG**