



DENTIRO® Wipes

Fully virucidal disinfectant wipes with a low alcohol content for small surfaces.

- **Kills 101 microorganisms**
- **Active in 60 seconds**
- **Compatible with more than 275 materials including artificial leather**
- **Residue-free drying**

DENTIRO® Wipes are ready-to-use disinfectant wipes for the residue-free disinfection and cleaning of small surfaces of non-invasive medical devices such as operating tables, gurneys, IV poles, and dental chairs. The highly saturated wipes have a low alcohol content of less than 50%, which not only reduces the potential for allergies but also increases material compatibility. DENTIRO® Wipes are aldehyde-free and biodegradable.

Scope of Application



- For the fast-acting wipe disinfection and cleaning of small surfaces of medical devices, e.g. operating tables, suction hoses, handpieces, contra-angles, operating light handles, and medical equipment made of aluminum, copper, silicon, rubber, and plastics such as PVC and acrylic glass

Use



Completely moisten the surface or object to be disinfected and leave to act. Adhere to exposure times. If necessary, dry the surface or object with a clean single-use cloth. Use on alcohol-compatible materials only. Test before use on sensitive materials by applying to a small inconspicuous area.



Effective Spectrum / Exposure Times



Bactericidal	15 seconds
Antibiotic-resistant bacteria	15 seconds
Tuberculocidal	60 seconds
Mycobactericidal	60 seconds
Yeasticidal	15 seconds
Antibiotic-resistant yeasts	15 seconds
Enveloped viruses	15 seconds
BVDV	15 seconds
Coronavirus	15 seconds
Hepatitis B virus	15 seconds
Hepatitis C virus	15 seconds
Herpes simplex virus	15 seconds
HIV	15 seconds
Influenza A virus	15 seconds
Vaccinia virus	15 seconds
Non-enveloped viruses	60 seconds
Adenovirus	60 seconds
Norovirus	60 seconds
Poliovirus	60 seconds
Rotavirus	30 seconds

Complies with EN 13624, EN 13727, EN 14348, EN 14476, EN 16615 and RKI/DVV guidelines. Dermatologically tested.





Compatibility



Polyvinyl chloride	Compatible.
Silicone rubber	Compatible.
Rubber	Compatible.
Aluminium	Compatible.
Copper	Compatible.
Brass	Compatible.
Acrylic glass	Compatible.
Polyethylene	Compatible.
Polytetrafluoroethylene	Compatible.
Polypropylene	Compatible.
Cast-iron	Compatible.

Warning



H228 Flammable solid. H319 Causes serious eye irritation. P210 Keep away from heat/sparks/open flames/hot surfaces. No smoking. P233 Keep container tightly closed. P280 Wear protective gloves and eye protection. P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. P337+P313 If eye irritation persists: Get medical advice/attention. P501 Dispose of contents/container in accordance with relevant regulations.

Composition of the moisturizing solution in 100 g



31.3 g ethanol, 16.7 g isopropanol

Physical Properties



Physical state:	Liquid absorbed on fleece
Colour:	Colourless
Scent:	Of alcohol
pH-value:	Neutral
Shelf Life:	36 months (unopened container)

Delivery Forms



Refill pouch with 120 disinfectant wipes (14.5 x 20 cm), classic scent, Refill pouch with 120 disinfectant wipes (14.5 x 20 cm), lemon scent, Refill pouch with 120 disinfectant wipes (14.5 x 20 cm), floral scent, DENTIRO Wipes Tub T1 (Empty tub), (CH-041079, CH-041083, CH-041084, CH-042081)

Your Local Distributor



OCC Switzerland

Oro Clean Chemie AG
Allmendstrasse 21
8320 Fehraltorf
Switzerland
Telephone: +41 (0)44 226 44 44
Fax: +41 (0)44 226 44 00
info@oroclean.com
www.oroclean.com
<https://www.facebook.com/Oroclean>