

Operating Instructions

Status of technical operating instructions: 01/96

1. General terms:

- 1.1. Product description:** Sati Brush No. 310-0000
- 1.2. Delivery form:** Polishing wheels in a small plastic bag Content 50 pc.
diameter 13 mm
- 1.3. Manufacturer:** YETI Dentalprodukte GmbH - Industriestraße 3 - D-78234 Engen/Germany
Tel.: 0049-7733/9410-0 FAX: 0049/7733-9410-22
- 1.4. Storage:** - store dry
- protect from blow or other impact

2. Use according to determinations:

Polishing wheels for treating of metal surface outside of patients mouth. Metal oxide can be removed from metal surface.

For a mat finish in order to recognize occlusal contact spots in the phase of provisional wearing of denture.

Do not use in patients mouth !

Any claims on warranty or liability will not be accepted if not used according to determinations !

3. Risk and danger during application:

1. Screw the Sati Brush wheel onto suitable mandrel only.
2. Use on rotation speed of 5.000-10.000 per minute.
3. Keep away from unqualified persons (e.g. children) because of the danger of swallowing. Only to be used by authorised persons.
4. Use only with safety glasses and a mouth protection.

4. Operating Instructions:

4.1. Preparation

Screw the Sati Brush wheel onto a suitable mandrel according to the speed of rotation.

4.2. Application

Satine (mat) the desired surface with moderate pressure observing the suggested engine speed range.

Hold the object which is to be satined firmly with your hand.

The polishing effect is less when Sati Brush wheel is worn-out. Change Sati Brush wheel when needed.

4.3. After application:

- Clean the polished object dry or with water.
- Sati Brush can be used again.

5. Disposal: According to current local regulations.

6. Durability: Consumer goods should be used rapid after opening and used up within a few months.

