

Operating Instructions

Status of technical instructions: 05/01

1. General terms:

- 1.1. Product description:** **DIE SPACERS (CLEAR SPACER/ COLOR SPACER)**
several colours/ thicknesses
THINNER for spacers

1.2. Delivery includes:

CLEAR SPACER transparent, 3my	#501-0001	1 bottle/	#501-0006	6 bottles of 20ml
COLOR SPACER red transparent, 7my	#500-0001	1 bottle/	#500-0006	6 bottles of 20ml
DIE SPACER yellow, 7 my	#502-0001	1 bottle/	#502-0006	6 bottles of 20ml
DIE SPACER blue, 10 my	#505-0001	1 bottle/	#505-0006	6 bottles of 20ml
DIE SPACER dentin colour, 10 my	#506-0001	1 bottle/	#506-0006	6 bottles of 20ml
DIE SPACER gold, 13 my	#504-0001	1 bottle/	#504-0006	6 bottles of 20ml
DIE SPACER silver, 13 my	#503-0001	1 bottle/	#503-0006	6 bottles of 20ml
THINNER for spacers	#510-0001	1 bottle/	#510-0006	6 bottles of 20ml

- 1.3. Manufacturer's address:** YETI Dentalprodukte GmbH - Industriestrasse 3 - D-78234 Engen
phone: 0049-7733-9410-0 FAX: 0049-7733-941022

- 1.4. Storage conditions:** Store in dry and cool place. Keep container closed.

2. Application according to determinations:

YETI SPACERS are to be applied on plaster dies prior to the wax sculpturing. SPACERS seal the plaster and serve to simulate the space needed for fixing cement between the crown and the die. The quality of YETI SPACERS enables robust, even and homogenous layers.

Any claims on warranty or liability will not be accepted if not used according to determinations!

3. Caution:

Ensure that unqualified persons (e.g. children) do not come into contact with the product.. For dental technicians and authorised persons only. Close bottles after use in order to avoid inhalation or drying out of the spacers.

4. Operating instructions:

4.1. Preparation:

Clear and red Color Spacer: No need to shake the bottle./ Other colours: Shake the bottle well before use. Open the bottle and take out the brush carefully.

4.2. Application:

Clear Spacer: Apply onto the entire plaster die including the area below the margin line. A layer of 3 my on dry plaster will seal the surface of the die as a base for a chemical bonding with further spacers. Other spacers: Apply the chosen die spacer (according to color/ thickness, see point 1.2) 1 mm above the margin line. After drying, several layers can be applied to achieve the ideal thickness.

4.3. After use:

After application the bottle must be closed immediately.

5. Disposal:

According to current local regulations.

6. Durability:

Material should be used within few months after opening.