Description

Multilink Automix try-in pastes are water-soluble glycerine pastes. The try-in pastes enable the simulation of the shade effect of all-ceramic restorations (inlays, onlays, crowns) that are to be cemented with Multilink Automix luting composite.

It is important to note that the shade and translucency of the try-in pastes corresponds to that of the **cured (!)** Multilink Automix luting composite.

Shades

Multilink Automix try-in pastes are available in four shades with different degrees of translucency:

- transparent (high translucency)
- yellow (high translucency)
- opaque (low translucency)
- white (medium translucency)

The shade of the try-in paste may subtly differ from that of the corresponding uncured Multilink Automix luting composite, since the luting composite attains its final shade only after it has been cured.

The shade effect of the try-in pastes is coordinated with that of the fully cured luting composite of the the same thickness.

Composition

The Multilink Automix try-in pastes contain glycerine, mineral fillers and dyes.

Indication

Simulation of the shade of the permanent luting composite.

Contraindication

Do not use if the patient is known to be allergic to any of the ingredients of the Multilink Automix try-in pastes.

Side effects

The hygroscopic try-in pastes cause fluid movement within the dentin tubules. As a result, the patient may feel temporary discomfort. Therefore, the try-in procedure should be carried out under anaesthetic or after desensitization with a desensitizer (e.g. Telio Desensitizer).

Interactions

Completely rinse any try-in pastes residue from the restoration and the tooth with water prior to applying the permanent cement. Otherwise, the bond between the restoration and Multilink Automix may be compromised.

Application

- Shade simulation with the help of try-in pastes must always be carried out before the working field is isolated or the teeth are dried, since drying temporarily lightens the tooth structure.
- Apply a small amount of the selected try-in paste to the restoration and carefully seat it without pressure. If the restoration cannot be seated easily, check the inner surface of the restoration as well as the proximal regions for any premature contacts. Do not force the restoration into place (risk of fracture!).

Note:

Care should be taken when checking the occlusion of brittle ceramic restorations before they are permanently cemented, as there is a risk of fracture.

- The try-in pastes are water-soluble and can be easily removed under running water or with ultrasonic cleaners. The restoration should be thoroughly cleaned and dried.
- Thoroughly clean the tooth structure with water spray.
- Permanently seat the restoration with Multilink Automix (See Multilink Automix Instructions for Use).

Storage

- Do not use Multilink Automix try-in pastes after the expiry date.
- Store at 2-28 °C / 36-82 °F.
- Expiry date: See information on syringes and packaging

Keep material out of children's reach.**For use in dentistry only.**

This material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damage resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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www.ivoclarvivadent.com

Ivoclar Vivadent Pty. Ltd.

1 – 5 Overseas Drive, P.O. Box 367, Noble Park, Vic. 3174, Australia

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www.ivoclarvivadent.com

Ivoclar Vivadent Marketing Ltd.

Rm 603 Kuen Yang, International Business Plaza, No. 798 Zhao Jia Bang Road,

Shanghai 200030, China

Tel. +86 21 5456 0776, Fax +86 21 6445 1561

www.ivoclarvivadent.com

Ivoclar Vivadent Marketing Ltd.

Calle 134 No. 7-B-83, Of. 520, Bogotá, Colombia

Tel. +57 1 627 33 99, Fax +57 1 633 16 63

www.ivoclarvivadent.com

Ivoclar Vivadent SAS

B.P. 118, F-74410 Saint-Jorioz, France

Tel. +33 450 88 64 00, Fax +33 450 68 91 52

www.ivoclarvivadent.fr

Ivoclar Vivadent GmbH

Dr. Adolf-Schneider-Str. 2, D-73479 Ellwangen, Jagst, Germany

Tel. +49 (0) 79 61 / 8 89-0, Fax +49 (0) 79 61 / 63 26

www.ivoclarvivadent.de

Ivoclar Vivadent Marketing (India) Pvt. Ltd.

503/504 Raheja Plaza, 15 B Shah Industrial Estate, Veera Desai Road, Andheri (West),

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www.ivoclarvivadent.com.mx

Ivoclar Vivadent Ltd.

12 Omega St, Albany, PO Box 5243 Wellesley St, Auckland, New Zealand

Tel. +64 9 914 9999, Fax +64 9 814 9990

www.ivoclarvivadent.co.nz

Ivoclar Vivadent Polska Sp. z o.o.

Al. Jana Pawla II 78, PL-00175 Warszawa, Poland

Tel. +48 22 635 54 96, Fax +48 22 635 54 69

www.ivoclarvivadent.pl

Ivoclar Vivadent Marketing Ltd.

Derbenevskaja Naberezhnaya 11, Geb. W, 115114 Moscow, Russia

Tel. +7 495 913 66 19, Fax +7 495 913 66 15

www.ivoclarvivadent.ru

Ivoclar Vivadent Marketing Ltd.

Qlaya Main St., Siricon Building No.14, 2nd Floor, Office No. 204,

P.O. Box 300146, Riyadh 11372, Saudi Arabia

Tel. +966 1 293 83 45, Fax +966 1 293 83 44

www.ivoclarvivadent.com

Ivoclar Vivadent Marketing Ltd.

171 Chin Swee Road, #02-01 San Centre, Singapore 169877

Tel. +65 6535 6775, Fax +65 6535 4991

www.ivoclarvivadent.com

Ivoclar Vivadent S.L.U.

c/ Emilio Muñoz Nº 15, Entrada c/ Albarracin, E-28037 Madrid, Spain

Tel. + 34 91 375 78 20, Fax + 34 91 375 78 38

www.ivoclarvivadent.es

Ivoclar Vivadent AB

Dalvägen 14, S-169 56 Solna, Sweden

Tel. +46 (0) 8 514 93 930, Fax +46 (0) 8 514 93 940

www.ivoclarvivadent.se

Ivoclar Vivadent Liaison Office

Ahi Evran Caddesi No 1, Polaris Is Merkezi Kat: 7, 80670 Maslak, Istanbul, Turkey, Tel.

+90 212 346 04 04, Fax +90 212 346 04 24

www.ivoclarvivadent.com

Ivoclar Vivadent Limited

Ground Floor Compass Building, Feldspar Close, Warrens Business Park,

Enderby, Leicester LE19 4SE, United Kingdom

Tel. +44 116 284 78 80, Fax +44 116 284 78 81

www.ivoclarvivadent.com

Ivoclar Vivadent, Inc.

175 Pineview Drive, Amherst, N.Y. 14228, USA

Tel. +1 800 533 6825, Fax +1 716 691 2285

www.ivoclarvivadent.com