
	TECHNICAL DATA SHEET	CODE: TDS_095
	REFRIGERATING SPRAY PROCLINIC	DATE: 26/JUN/19
	PROCLINIC	VERSION: 01

PRODUCT	REFRIGERATING SPRAY PROCLINIC	
REFERENCE	01924 01923	
COMPANY	PRISMAN GmbH. Otto-Hahn-Ring 6-18 64653 Lorsch Tel: +49 6251 866980-0	
CLASSIFICATION	CLASS I MEDICAL DEVICE Directive 93/42/CEE	
PRODUCT DESCRIPTION	Available in two pleasant fragrances such as orange and mint, it is applied with precision thanks to its special applicator that also saves spray in each application. Cool up to - 40 ° C.	
SPECIFICATIONS	01924: Mint Fragrance 01923: Fragrance Orange	
PACKAGING	Bottle of 200 mL.	
SHELF LIFE	3 years	
INTENDED OF USE	Coolant spray for pulp vitality test and cooling of pellets, glass tiles, etc.	
SECURITY	Keep away from heat, hot surfaces, sparks, flames and other sources of ignition. Do not smoke nearby. Do not pierce or burn, even after use. Do not spray on an open flame or other source of ignition. Protect from sunlight. Do not expose to temperatures above 50 ° C.	
TRANSPORT AND STORAGE CONDITIONS	Store in a cool place. Keep container tightly closed. Store in a cool, dry place in tightly closed containers. Protect from heat and direct sunlight. Store in vertical position.	