

COLORED PAPER/PLASTIC AND PLASTIC DENTAL APRONS IN ROLL

Disposable protective aprons in rolls are available in two qualities: paper/plastic and plastic.

The aprons in roll are made by 1-2 plies of pure long fiber cellulose wadding coupled with one ply of impermeable high density polyethylene. Plastic aprons in roll are made by impermeable high density polyethylene.

The product perfect absorbency and water resistance ensures a high patient comfort.

Printing of own brand is possible with minimum order to discuss. Lot number is printed on outer carton.

They are available with the following sizes and colours:

1 ply paper + 1 ply polyethylene		
MA010160	80 sheets cm.48x60 light blue	12 aprons / carton
MA010260	80 sheets cm.48x60 green	12 aprons / carton
MA010360	80 sheets cm.48x60 yellow	12 aprons / carton
MA010460	80 sheets cm.48x60 blue	12 aprons / carton
MA010560	80 sheets cm.48x60 lavender	12 aprons / carton
MA010660	80 sheets cm.48x60 peach	12 aprons / carton
MA011060	80 sheets cm.48x60 pink	12 aprons / carton
MA010860	80 sheets cm.48x60 lime	12 aprons / carton
MA010960	80 sheets cm.48x60 orange	12 aprons / carton
MA011160	80 sheets cm.48x60 white	12 aprons / carton
1 ply paper + 1 ply polyethylene		
MA010180	60 sheets cm.48x80 light blue	12 aprons / carton
MA010280	60 sheets cm.48x80 green	12 aprons / carton
MA010380	60 sheets cm.48x80 yellow	12 aprons / carton
MA010480	60 sheets cm.48x80 blue	12 aprons / carton
MA010580	60 sheets cm.48x80 lavender	12 aprons / carton
MA010680	60 sheets cm.48x80 peach	12 aprons / carton

MA011080	60 sheets cm.48x80 pink	12 aprons / carton
MA010880	60 sheets cm.48x80 lime	12 aprons / carton
MA010980	60 sheets cm.48x80 orange	12 aprons / carton
MA011180	60 sheets cm.48x80 white	12 aprons / carton
<u>Plastic</u>		
MA020160	200 sheets cm.48x60 blue	12 aprons / carton
MA020260	200 sheets cm.48x60 green	12 aprons / carton
MA020360	200 sheets cm.48x60 white	12 aprons / carton
<u>Plastic</u>		
MA020180	150 sheets cm.48x80 blue	12 aprons / carton
MA020280	150 sheets cm.48x80 green	12 aprons / carton
MA020380	150 sheets cm.48x80 white	12 aprons / carton

STORAGE CONDITIONS: keep at room temperature in a cool and dry place, far from sunlight and fire.

Ospitaletto (Brescia) 26.02.2016

Edition 1

SOGEVA S.r.l.
Via I°Maggio traversa terza n° 18/22
25036 OSPITALETTO (Brescia)
Tel. 030.2500251 - Fax 030.2500482
Codice Fiscale 09056160154
Partita I.V.A. 03558360172