
	TECHNICAL DATA SHEET	CODE: BD26801
	PROPHY CUP	DATE: 10/Mar/19
	BESTDENT	VERSION: 01

<b>PRODUCT</b>	Prophy Cup	
<b>REFERENCE</b>	BD26801	
<b>COMPANY</b>	SHANGHAI DOCHEM INDUSTRIES CO., LTD.	
<b>CLASSIFICATION</b>	CLASS IIa MEDICAL DEVICE Directive 93/42/CEE	
<b>PRODUCT DESCRIPTION</b>	Latex Free Smooth soft, with high torque gear to perform, will not skip or over heat, long lasting.	
<b>SPECIFICATIONS</b>	Latch type, 4 webs	
<b>PACKAGING</b>	Packing: 100pcs/box, 100boxes/ctn	
<b>SHELF LIFE</b>	3 years	
<b>INTENDED OF USE</b>	Used for polishing and removal of heavy stains	
<b>SECURITY</b>	We finish the risk analysis according to standard ISO14971. According to the identification of the potential risk and the assessment of the every risk, then use the still of FMEA to evaluate the R.P.N., the no. is very low. It means, the product has high security. To user, it is an acceptable risk causing by every potential danger.	
<b>TRANSPORT AND STORAGE CONDITIONS</b>	The transport and storage is in accordance with the relevant requirements after assessment, the mode of transport and storage is reasonable and feasible, can be used for the usual production process totally.	



Signature/Stamp:

Date: March 10th, 2020