
	TECHNICAL DATA SHEET	CODE: TDS_037
	COMPOSITE PM MICROHIBRIDO FOTOPOLIMERIZABLE	DATE: 26/JUN/19
	BESTDENT	VERSION: 01

PRODUCT	MICRO- HYBRID LIGHT CURE COMPOSITE A1, A2, A3, A3.5 (600-77, 78, 79, 80) 4g Syringe	
REFERENCE	60077 60078 60079 60080	
COMPANY	CONAMCO José María Bustillos 20, Algarín, 06880 Ciudad de México, México	
CLASSIFICATION	CLASS IIa MEDICAL DEVICE Directive 93/42/CEE	
PRODUCT DESCRIPTION	Physical characteristics: Soft paste, Sensitive to light. Its mechanical properties as well as its brightness and wear resistance is very similar to those found in the market.	
PRODUCT COMPOSITION	Composition: Crystalline and smoked silicon dioxide, methacrylates, stabilizers and photo-initiators. Microhybrid Composite - This composite has the correct proportion of 77% microparticles, all of them smaller than 5 microns.	
SPECIFICATIONS	Photopolimerizable composite, only for dental professionals.	
PACKAGING	1 Syringe of 4g	
SHELF LIFE	3 years	
INTENDED OF USE	Microhybrid composite for aesthetic restorations in previous ones. It also allows restorations in molars when aesthetic demands require it.	
SECURITY	Slight irritating to eyes and skin. May cause Sensitization by skin contact. Only adequately trained personnel should handle this product. This product is intended for use in dentistry only.	
TRANSPORT AND STORAGE CONDITIONS	Store only in the original package. Keep package tightly sealed. Store in a dry area. Store away from food and beverages. Protect from excessive heat and cold and direct sunlight.	

N/A: Not Aplicable