



Com-Fit™ Masks

Instructions For Use



1. Perform hand hygiene before touching mask.



2. Orient mask for proper placement, pleats down.



3. Holding mask by earloops, place loops around each ear.



4. Mold malleable strip at the top edge of mask to shape of nose.



5. Pull bottom of mask over mouth and chin.



Application Area

Com-Fit Masks are a complete line of ASTM-rated face masks, designed to provide maximum comfort and protection to meet the needs of any office.

ASTM-1

USAGE

For procedures with light fluid, aerosols and/or sprays.

PROCEDURES:

Patient Exams, Operatory Cleaning/ Maintenance, Impressions, Lab Trimming, Finishing & Polishing, Orthodontics.

ASTM-2

USAGE

For procedures with low to moderate fluid, aerosols and/or sprays.

PROCEDURES:

Restorative/Composites, Endodontics, Prophylaxis, Sealants, Scaling & Root Planing, Limited Oral Surgery.

ASTM-3

USAGE

For procedures with moderate to heavy fluid, aerosols and/or sprays.

PROCEDURES:

Crown Preparation, Implant Placement, Use of Ultrasonic Scaler, Use of Piezo Scaler with Water, Periodontal Surgery, Complex Oral Surgery.

Proper Fit + Proper Wearing = Maximum Protection



Com-Fit™ Masks

Filtration, breathability and comfort.



Better

- More breathable than leading competitor masks.¹
- More comfortable than leading competitor masks.¹



Safer

- Better fit than leading competitors for improved protection.¹
- Exceeds ASTM performance requirements for bacterial and particle filtration.¹
- Patented nose bend for proper mask placement.



Faster

- Only mask dispenser with large side panel opening for quick and easy access to masks.



Com-Fit Plush
Natural Fit Masks
ASTM 1



Com-Fit Plush
Natural Fit Masks
ASTM 2



Com-Fit Plush
Natural Fit Masks
ASTM 3

Packaging & Ordering Information

20355	Plush Natural Fit ASTM 1, box of 40
20356	Plush Natural Fit ASTM 2, box of 40
20357	Plush Natural Fit ASTM 3, box of 40
20352	Plush Masks ASTM 1, box of 40
20353	Plush Masks ASTM 2, box of 40
20354	Plush Masks ASTM 3, box of 40
20344	Super Sensitive Tie On White (die free), box of 50

20346	Super Sensitive Ear Loop White (die free), box of 50
20326	Super High Filtration Blue Ear loop, box of 40
20336	Super High Filtration Tie On Blue, box of 40
20350	Groovy Mask, box of 40
20324	Fluid Resistant Proc Blue, box of 50
20325	Fluid Resistant Proc Pink, box of 50
20340	Easy Breathe Ear Loop Blue, box of 40
20342	Easy Breathe Tie On, box of 40

For more information: dentsplysirona.com/preventive.

© 2018 Dentsply Sirona Preventive. 0010322LI Rev. 1 (0118)

¹ Data on File