



TECHNICAL DATA

Brand name: Hoffmann's Impression Compound green

Generic name: Impression compound

Purpose: Bite registration; copper-ring impressions; correction of models; raising of impression edges.

Use Form and Packaging (size, size of all variations):

Package with 15 sticks, 135 g

Composition:

stearine
copal
talcum
carnauba wax
colour

Toxicological data:

Contra-indication:	None
Side effects:	None
Precautions:	None

According to our experiences there are no noxious effects known when used and handled according to specifications.

Safety data sheets and instructions available at www.hoffmann-dental.com.



Hoffmann's Impression Compound green
Technical data page 2/2

Other data:

Detailed data on the submitted certificates
Certificate DIN EN ISO 13485:2016
Certificate Annex II of the council directive 93/42/EEC concerning medical devices

UMDNS-No.: 16-680

Customs Tariff Number:
3407 0000

Shelf Life: 60 month


Storage conditions: 2–28 °C, dry

Catalogue number:
82607

Date:

10.05.2020

Signature:



Yvonne Hoffmann
(Managing Director)