



## 5. i Led 产品规格

### i Led Product Specification

#### 1、机器正常工作条件:

Working condition:

1.1、温度:  $5^{\circ}\text{C} \sim 40^{\circ}\text{C}$ 。

1.1 Working temperature:  $5^{\circ}\text{C} \sim 40^{\circ}\text{C}$ 。

1.2、相对湿度: 10%~93%。

1.2 Relative humidity: 10%~93%。

#### 2、贮存条件:

2.1、产品应贮存在于相对湿度不超过 10%~93%,大气压力为  
 $70\text{KPa} \sim 106\text{KPa}$ ,  $-20^{\circ}\text{C} \sim +55^{\circ}\text{C}$ ,通风良好的室内。

Storage conditions:

Relative humidity: 10%~93%

Atmosphere pressure:  $70\text{kPa} \sim 106\text{kPa}$

Environment temperature:  $-20^{\circ}\text{C} \sim +55^{\circ}\text{C}$

2.2、贮存时不与有毒、有腐蚀性、易燃、易爆的物品混放

Don't put it together with poisonous、inflammable、explosive、substances  
and corrosive materials during transportation

#### 3、机器各参数:

Main technical specifications



### 3.1、电池：可充电锂离子电池

Battery: rechargeable Lithium battery

电池型号： ICR 18490

Battery model: ICR 18490

标称电压： 3.6V

Battery voltage: 3.6V

标称容量： 1400mAh

Battery capacity: 1400mAh

### 3.2、消耗功率：不大于 8W

Consumption power: ≤8W

### 3.3、光源：（1）蓝光 （2）波长：385nm-515nm。

Light source: (1) blue light (2) wave length: 385nm-515nm。

### 3.4、外形尺寸：240mm×25mm×25mm(长×宽×高)

Dimensions: 240mm×25mm×25mm

### 3.5、净重：278g

Net weight: 278g

### 3.6、充电器工作电压：AC100—240V 频率：50/60Hz

Input of adapter: 100V- 240V AC 50/60HZ

### 3.7、充电器内置熔断器：T1A250V

Built-in fuse of adapter: T1A250V

File No.: ZMN/WI-07-025-05

Version: B

3.8、按防电击类型分类: II 类设备。

Type of protection against electric shock: Class II

3.9、按防电击的程度分类:B 型应用部分

Degree of protection against electric shock: Type B applied part

3.10、按对进液的防护程度分类:普通器材(IPX0)。

Degree of protection against harmful ingress of water: Ordinary equipment (IPX0)

3.11、易燃麻醉混合物和空气,氧气或一氧化二氮的安全度数:

不适合在易燃麻醉混合物和空气,氧气或一氧化二氮中使用的器材。

Degree of safety of application in the presence of a flammable anesthetic being used in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.