

Core build-up Proclinic Expert



Core build-up **Proclinic Expert** is a dual-cured composite for performing core build-ups, which is highly radiopaque and flowable with excellent mechanical properties. Presented in an automix syringe for easy and convenient application.

COMPRESSIVE STRENGTH

Its high compressive strength allows the creation of very durable restorations.

SC: self-curing, LC: light-curing*

Core build-up Proclinic Expert

322 MPa (SC) → 349 MPa (LC)

Other brands

250 MPa → 300 MPa

FILLER LOADING

The high filler loading of Proclinic Expert's core build-up, 70% by weight, ensures high wear resistance.

DENTINE BONDING STRENGTH

High bonding strength to dentine (22 Mpa), which ensures long-lasting restorations.

- Two indications in one product: core build-ups and post cementation.
- Good thixotropic properties ensuring good product handling and easy adaptation to irregular preparations and endodontic posts.
- High radiopacity which allows long-term treatment monitoring.
- Efficient mixing tips which deliver exactly the required amount.
- Optimal mixing tips to facilitate posterior access.

Proclinic Expert's core build-up
1 x 10g syringe, 5 mixing tips and 5
intra-oral tips (type 4)

REORDER
78560





DENTINE SHADE

The dentine shade allows you to create highly aesthetic, non-metal restorations.

SOLUBILITY IN WATER

Achieve longer lasting restorations thanks to its low water solubility.

Core build-up Proclinic Expert

10 $\mu\text{g}/\text{mm}^3$

Other brands

25 $\mu\text{g}/\text{mm}^3$

WATER ABSORPTION

Its low water absorption helps to create stable, long-term restorations with no risk of discolouration.

Core build-up Proclinic Expert

0.1 $\mu\text{g}/\text{mm}^3$

Other brands

2.5 $\mu\text{g}/\text{mm}^3$



Dr. Eduardo González Arroyo

Clínica Fernández

www.oralsturgerytube.com

Member of Association of Dental Surgeons number 28007247


"The intrinsic problems faced when doing core build-ups are overcome by using a strong, easy-to use and fast setting material, such as the new product Proclinic Expert's core build-up"

Dr. Carlos Gargallo

Clínica Carlos Gargallo

Member of Association of Dental Surgeons number 3947

"Using Proclinic Expert's core build-up is very satisfying, as it has a very good viscosity, withstands preparation and is easy to apply. Furthermore, the change of colour on setting indicates correct polymerisation."



Proclinic Expert's core build-up

REORDER
78560

