

RESTORATIONS



Compressive
strength
400MPa

61% filler
loading

59.3% opacity

2.05%
volumetric
shrinkage

SIMPLE AND LONG-LASTING RESTORATION

Rapid filling,
minimal shrinkage stress,
optimal precision

ANTERIOR AND POSTERIOR



Aesthetic precision in **2** mm increments

Proclinic Expert composites are a perfect material
for both anterior and posterior restorations

Dr Lotte Hallander

TandCity ApS, København, Denmark.
www.tandcity.dk

"Handling and manageability of Proclinic Expert Nano Hybrid is very easy and the nano effect provides very good shade matching. I truly recommend this product."



Dr Alexander Johansson

Clinique ZahnCity GmbH, Zurich,
Switzerland.
www.zahncity.ch

"The Proclinic Expert Nano Hybrid composite has optimal viscosity and is easily moulded and shaped. It is very flexible, wear-resistant and has excellent physical strength."



Rapid filling in layers of up to **5** mm thick

The **Proclinic Expert** Bulk Fill composite
is suitable for treating large cavities
in posterior restorations

Opinions of several
consultants from the US-based
dental test institute,
"The Dental Advisor"



"It is ideal for rapid restorations in
posterior teeth without overlays"

"It has a high curing depth and
high strength"

"Great for contouring, finishing
and polishing"

Proclinic Expert Bulk Fill

The fastest restoration



Proclinic Expert Bulk Fill is a light-cured, radiopaque posterior composite capable of treating a maximum thickness of 5mm in 20 seconds.



- 5mm increments can be light-cured in 20 seconds.
- Does not stick to instruments when they are removed.

- Good viscosity, condenses well in posterior restorations.
- It is firm and can be condensed to create better contact points.

Comments from consultants
at the institute "The Dental Advisor"



"Proclinic Expert Bulk Fill has very good handling properties"

"Proclinic Expert Bulk Fill has a treatment depth of 5mm"

"Shading for posterior teeth is nonessential and Proclinic Expert Bulk Fill composite is very suitable for particularly large cavities"

Proclinic Expert Bulk Fill

Excelente nota
en los tests clínicos




Proclinic Expert Bulk Fill
+ + + + ½

THE DENTAL ADVISOR
www.dentaladvisor.com

RATINGS:
Excellent + + + + +
Very Good + + + +
Good + + +


345
CLINICAL
APPLICATIONS



91%

SATISFACTION

78%
of consultants consider
Proclinic Expert Bulk Fill
TO BE AS GOOD
OR BETTER
than other brands



Proclinic Expert Nano Hybrid

Your ideal anterior/posterior composite




Proclinic Expert Nano Hybrid composite is perfect for both anterior and posterior restorations. It contains a fluorescent pigment that ensures the natural appearance of teeth and is very easy to handle.

COLOUR STABILITY

Proclinic Expert Nano Hybrid is virtually unaffected by accelerated ageing.


*Change of colour from shade A2 after ageing in transmission mode**



TRANSLUCENCY

Thanks to its translucency, **Proclinic Expert Nano Hybrid** reflects the surrounding tooth structure, guaranteeing an invisible restoration.


*Opacity**



COMPRESSIVE STRENGTH

The high compressive strength of Proclinic Expert Nano Hybrid helps it withstand chewing forces, increasing its life expectancy.

*Strength**



* Results of testing by the manufacturer comparing the product with the three most popular brands available on the Spanish market.

Proclinic Expert Nano Flow

Your most versatile
composite



Proclinic Expert Nano Flow is indicated for almost all types of restorations. It is radiopaque and guarantees the natural appearance of teeth. It contains nanofillers that improve polishing properties, wear resistance and compressive strength, and UDMA resins that reduce viscosity and minimise post-operative sensitivity.

COMPRESSIVE STRENGTH

The high compressive strength of **Proclinic Expert Nano Flow** means it has a wider range of application possibilities.


MPa at 1 hour*



FILLER LOADING

The high filler loading of **Proclinic Expert Nano Flow**, 64% by weight, lends it an excellent wear resistance.

FLUORIDE RELEASE



Proclinic Expert Micro Hybrid

A light-cured, radiopaque, microhybrid composite for anterior/posterior restorations. It does not stick to instruments and has low shrinkage, which improves application and clinical results.

