

FICHA TÉCNICA

MECHERO DE ALCOHOL METÁLICO		
REFERENCIA	80221	
FABRICANTE / RESPONSABLE	MARK ONE SURGICAL Ugoki Road, Adalat Garah, Sialkot 51330, Pakistan	
DESCRIPCIÓN	Mechero de alcohol metálico con mecha incluida.	
ESPECIFICACIÓN	Diámetro: 7cm. Alto: 4 cm.	
COMPOSICIÓN	Acero inoxidable.	
ENVASE	1 unidad.	
FINALIDAD PREVISTA	- Se utiliza tanto en laboratorios como en clínicas dentales para realizar experimentos como calentar, esterilizar instrumentos, etc. - Principalmente es usado para calentar la cera y poderla manipular en la realización de prótesis dentales.	
ALMACENAMIENTO Y TRANSPORTE	No requiere condiciones especiales de almacenamiento y transporte.	

TECHNICAL DATA SHEET

ALCOHOL BURNER METALLIC		
REFERENCE	80221	
MANUFACTURER	MARK ONE SURGICAL Ugoki Road, Adalat Garah, Sialkot 51330, Pakistan	
DESCRIPTION	Metallic alcohol lighter with wick included.	
SPECIFICATION	Diameter: 7 cm. Height: 4 cm.	
COMPOSITION	Stainless steel.	
PACKAGE	1 unit.	
INTENDED USE	- It is used both in laboratories and in dental clinics for experiments such as heating, sterilizing instruments, etc. - It is mainly used to heat the wax to be able to manipulate it in the production of dental prostheses.	
STORAGE AND TRANSPORT	It does not require special storage and transportation conditions.	