


**MESTRA®**

FICHA TECNICA

Reference:	Ref. 080087	Picture: 
Name:	PLASTIC DECANter WITH COLLECTION BAG SMALL SIZE	
Manufacturer:	TALLERES MESTRAITUA S.L. Txori-Erri Etorbidea, 60 Tel +34 944530388 Fax +34 944711725 info@mestra.es www.mestra.es 48150 SONDIKA – BILBAO SPAIN	
Description:	Made of resistant plastic. The design is specifically studied to achieve quick and easy cleaning because the dirt is collected in a plastic bag. The hermetic closure of its lid prevents the spread of bad odours. Its size and shape have been designed for easy installation under the sink.	
Replacements:	Plastic bags 20 units (ref. 080085-02)	
Technical data: Height: 345 mm Diameter: 350 mm Capacity: 14 L Weight: 2.2 kg		 The technical drawing includes a top view showing a circular lid with a diameter of 352 mm and a central hole with a diameter of 99.5 mm. The side view shows the decanter's height dimensions: a total height of 324 mm, a height to the top of the lid of 307 mm, and a height to the center of the lid of 241 mm.
Security measures:	• Choose a rigid, flat, and well-leveled surface for the location of the decanter. • Avoid the presence of electrical cables or circuits in the vicinity of the decanter. The minimum recommended safety distance should be at least 1 m. • Carry out the transport and movement operations of the decanter very carefully. Avoid risks of the decanter tipping over. • In case of difficulty in finding the optimal way to adapt the decanter to your laboratory's drainage system, consult a professional plumber.	