


**MESTRA®**

FICHA TECNICA

Reference:	Ref. 020100	
Bame:	TRIANGULAR FLASK	
Standards:	-	
Manufacturer:	TALLERES MESTRAITUA S.L. Txori-Erri Etorbidea, 60 Tel. + 34 944530388 Fax + 34 944711725 mestra@mestra.es www.mestra.es 48150 SONDIKA – BILBAO SPAIN	
Picture:		
Description:	<p>Made of brass.</p> <p>Consisting of four parts to facilitate removal of plaster block.</p> <p>Weight: 1.7 kg</p> <p>Inside dimensions:</p> <div data-bbox="740 1778 1026 2024"></div>	