


FICHA TÉCNICA

CERA PACIENTE ORTODONCIA PROCLINIC		
REFERENCIAS	L50933: Sin aroma L50934: Aroma menta	
FABRICANTE / RESPONSABLE	CHL MEDICAL SOLUTIONS, Srl Via Felice Cavallotti, 15 20122 Milano - Italia	
CLASIFICACIÓN	Producto Sanitario Clase I, según Reglamento (UE) 2017/745	
DESCRIPCIÓN	Barras de cera protectora para evitar el roce de los aparatos de ortodoncia	
COMPOSICIÓN	Mezcla de ceras y resinas /aceites parafínicos	
ESPECIFICACIONES	L50933 (Sin Aroma) L50934 (Aroma menta)	Color: Blanco Tamaño de las barras: 50 x 5 mm Peso: 4 g (estuche de 5 barras)
ENVASE	Caja de 36 estuches (5 barras por estuche)	
VIDA ÚTIL	10 años a partir de la fecha de fabricación	
FINALIDAD PREVISTA	Protección al roce de aparatos de ortodoncia.	
INFORMACIÓN SOBRE LA SEGURIDAD	Este producto no representa ningún peligro si se utiliza para su uso normal.	
ALMACENAMIENTO Y TRANSPORTE	Almacenar en lugar fresco y seco. Mantener alejado de la luz solar. Mantener alejado de toda fuente de ignición. Evitar los agentes oxidantes. No utilizar si el envase está dañado.	

TECHNICAL DATA SHEET

ORTHODONTIC PATIENT WAX PROCLINIC		
REFERENCES	L50933: Scent-free L50934: Mint -scented	
MANUFACTURER	CHL MEDICAL SOLUTIONS, Srl Via Felice Cavallotti, 15 20122 Milano - Italy	
CLASSIFICATION	Medical Device, Class I, according to Regulation (EU) 2017/745	
DESCRIPTION	Protective wax to prevent chafing from orthodontic appliances	
COMPOSITION	Mixture of waxes + resins/paraffinic oils.	
SPECIFICACIONES	L50933 (Scent-free) L50934 (Mint-scented)	Colour: White Bar size: 50 x 5 mm Weight: 4.0 g (case of 5 bars)
PACKAGING	Box with 36 cases (5 bars per case)	
SHELF LIFE	10 years from the date of manufacture	
INTENDED USE	It is a protective wax to prevent chafing from orthodontic appliances	
SAFETY INFORMATION	This product is not expected to present any unusual hazards in proper use.	
HANDLING AND STORAGE	Store in a cool dry place, away from direct sunlight. Keep away from heat and flames. Do not store near strong oxidizing agents. Do not use if the package is damaged.	