

	FICHA TÉCNICA	TDS_
	PUNTAS GUTTA PERCHA RECIPROC	Fecha: 01/JUN/18
	PROCLINIC	Versión: 01

PRODUCTO	PUNTAS GUTTA PERCHA RECIPROC	
----------	------------------------------	---

	TECHNICAL DATA SHEET	TDS_
	GUTTA PERCHA RECIPROC POINTS	Date: 01/Jun/18
	BESTDENT	Version: 01

PRODUCT	BESTDENT GUTTA PERCHA RECIPROC POINTS	DE Gutta-percha. Nur für den zahnärztlichen Gebrauch. In trockener und kalter Umgebung aufbewahren. Keinen Autoklav oder Hitze zum sterilisieren benutzen. EN Gutta-percha. For dental use only. Store in dry cool area. Do not use autoclave or heat sterilize. ES Gutta-percha. Exclusivamente para uso dental. Mantener en lugar fresco y seco. No utilizar el autoclave ni esterilizar con calor seco. FR Gutta-percha. Uniquement pour l'usage dentaire. Stocker dans endroit sec avec une température constante. Ne pas stériliser avec l'autoclave ou à travers d'une forte source de chaleur. IT Gutta-percha. Esclusivamente per uso dentale. Mantenere in luogo fresco e asciutto. Non utilizzare autoclave o in sterilizzatrici a caldo. PL Gutta-perka. Wyłącznie do użytku stomatologicznego. Przechowywać w suchym i chłodnym miejscu. Nie sterylizować w autoklawie ani nie poddawać sterylizacji termicznej. PT Gutta-percha. Apenas para uso dentário. Armazenar em área seca e fresca. Não usar esterilização em autoclave ou por calor seco. Contents: Gutta Percha, Zinc Oxide, Barium Sulphate, Colouring Agent. EC REP CHENNAI S.A. No. 15/16A Street No. 13, 8-10/15 Industrial Park, Tan Binh District, Hochiminh City - Vietnam Best Dent Co., Ltd. No. 15/16A Street No. 13, Industrial Park, Tan Binh District, Hochiminh City - Vietnam CE 0120 5
REFERENCES	01589 GUTTAPERCHA TAPER F1 PROCLINIC 01590 GUTTAPERCHA TAPER F2 PROCLINIC 01591 GUTTAPERCHA TAPER F3 PROCLINIC 01594 GUTTAPERCHA TAPER SURT F1-F2-F3 PROCLINIC 01595 GUTTAPERCHA TAPER SURT F4-F5 PROCLINIC	
COMPANY	Pearl Dent Co. Ltd. No3-15B, street No 13, Tan Binh Industrial Park, Hochiminh City, Vietnam.	
CLASSIFICATION	Class IIa medical device, Directive 93/42/EEC	
DESCRIPTION	Gutta Percha is a device made from coagulated sap of certain tropical trees intended to fill the root canal of a tooth. The Gutta Percha is softened by heat and inserted into the root canal, where it hardens as it cools.	
SPECIFICATIONS	Gutta Percha Point is round and have pink color, with a large head and a small head. Depending on the type and size, Gutta Percha Point has different sizes, which are classified according to each type as follows : - ISO Standard Size, Color coded #8 - #140 ; - Detrey Size (XXF, XF, F, M) ; - A, B, C, D Size ; - Taper 0.04/ Taper 0.06 ; - X-Prototype F1, F2, F3, F4, F5 ; - GC Type ; - Reguler Size ; - Japan Size ; - Microtip Size ; Proc size (R25, R40, R50)	
COMPOSITION	Zinc oxide, Barium sulfate, coloring agent	
PACKAGING	Packaging: Box of PS plastic: Slide box Quantification: Slide box (6 compartments) contain: 60 points / box; Inside of Slide box: 10 points / compartment.	
INTENDED USE	Only for dental use. Filling material (permanent use) of the cleaned root canals of tooth.	
TRANSPORT AND STORAGE CONDITIONS	Store in a cool dry area, away from direct sunlight. Protect from moisture. The device should not be sterilized by autoclave or heat.	