
	TECHNICAL DATA SHEET	
		DATE: 12/04/17

PRODUCT	ORAL TIPS	
REFERENCE	0124	
COMPANY/Manufacturer	SHANGHAI NSJ HARDWARE FACTORY Address: No.1118, Huateng Rd., Huaxin Town, Qingpu District, Shanghai, China, 201708.	
CLASSIFICATION	CLASS I MEDICAL DEVICE Directive 93/42/CEE	
PRODUCT DESCRIPTION	Made of food-grade PP For delivering the mixed dental impression material into the designated area.	
SPECIFICATIONS	4.2mm outlet ID	
PACKAGING	100pcs/bag	
SHELF LIFE	5 years	
INTENDED OF USE	Used together with yellow mixing tip, and delivering the mixed dental impression material into the designated area.	
SECURITY	Handling: No special measures necessary if handled according to instruction.	
TRANSPORT AND STORAGE CONDITIONS	Regular Transport, Normal temperature and pressure	

SHANGHAI NSJ HARDWARE FACTORY

Caroline Lee

April 12, 2017

