


PRODUCT	50m 1:1 Plastic Dispenser PROCLINIC	
REFERENCE	0156	
COMPANY	SHANGHAI NSJ HARDWARE FACTORY	
CLASSIFICATION	CLASS I MEDICAL DEVICE Directive 93/42/CEE	
PRODUCT DESCRIPTION	Item no.: CG-50-1-1 Plastic Dispenser Cartridge: two-component Capacity: For 50ml(1:1&2:1) Material: Reinforced Nylon & Metal Alloy Features & Benefits: Most durable, Save-effort, Excellent ergonomics, Accurate and consistent	
SPECIFICATIONS	50ml 1:1/2:1	
PACKAGING	Weight: 0.286kg / piece Package:50pcs/ctn Dimension: 610*420*400mm	
SHELF LIFE	5 years	
INTENDED OF USE	Used for dispensing 50ml 1:1 dental silicone impression materials, which is suitable for any-brand silicone.	
SECURITY	Yes	
TRANSPORT AND STORAGE CONDITIONS	Normal temperature and pressure.	