


REFERENCIA	16100232	CÓDIGO EAN	842965902069
PRODUCTO	Paño quirúrgico 21916 - (IS) 50x50cm		

COMPONENTES

<p>Paño quirúrgico 21916 - (IS) 50x50cm</p> <p>Paño desechable confeccionado en tela sin tejer compuesto de dos capas, que lo hacen impermeable y absorbente.</p> <p>Largo 50cm; Ancho 50cm.</p> <ul style="list-style-type: none"> - Tejido IS (56 g/m²): PP 30g laminado con PE 25μ. 	 <p>50cm</p> <p>50cm</p>
---	--

PRESENTACIÓN

Unidades embalaje 150

Presentación Estéril por óxido de Etileno

Alta 26/02/2013

Revisión 7/03/2018


NORMATIVA

Todos los productos de Cardiva Soluciones Integrales disponen de marcado CE, cumpliendo con la Directiva Europea 93/42/CE relativa a productos sanitarios.

Cardiva Soluciones Integrales tiene implantado y certificado un sistema de gestión de calidad y medio ambiente acorde a las normas ISO 9001; ISO 13485 e ISO 14001.

REFERENCE	16100232	BAR CODE	842965902069
PRODUCT	Surgical Drape 21916 - (IS) 50x50cm		

COMPONENTS

<p>Surgical Drape 21916 - (IS) 50x50cm</p> <p>Disposable drape is manufactured in nonwovens composed of two layers, one of which was waterproof fabric.</p> <p>Length 50cm; Width 50cm.</p> <ul style="list-style-type: none"> - IS Fabric (56 g/m²): laminated PP 30g, PE 25μ. 	 <p>50cm</p> <p>50cm</p>
--	--

PRESENTATION

Packing Units 150

Presentation Oxide Etylene

Entry date 26/02/2013

Review 7/03/2018

APPLICABLE STANDARD

All Cardiva Integral Solutions' products have CE mark according to European Directive 93/42/CE concerning Medical Devices.

Cardiva Integral Solutions has implemented and certified a quality and environment management system according to the standards ISO 9001; ISO 13485 and ISO 14001.