


FICHA TÉCNICA

MICROAPLICADORES BESTDENT		
REFERENCIAS	25650: Fino 25651: Regular	
FABRICANTE / RESPONSABLE	CHL MEDICAL SOLUTIONS, Srl Via Felice Cavallotti, 15 20122 Milán - Italia	
CLASIFICACIÓN	Producto Sanitario Clase I, conforme a la Directiva 93/42/CEE	
DESCRIPCIÓN	Los microaplicadores desechables son productos sanitarios invasivos de uso transitorio. Se utilizan para la aplicación de materiales en la cavidad oral. Tamaño fino (color amarillo) y regular (color verde).	
COMPOSICIÓN	Constituidos por unos microcepillos de fibras y mango de Polipropileno.	
EMBALAJE	Envase con 100 unidades	
VIDA ÚTIL	3 años a partir de la fecha de fabricación	
FINALIDAD PREVISTA	Los microaplicadores se utilizan para la aplicación de pequeñas cantidades de materiales (selladores, agentes adherentes, cementos...) en zonas de la cavidad bucal de acceso limitado.	
INFORMACIÓN SOBRE LA SEGURIDAD	Solamente para uso profesional dental. Producto de un solo uso. No estériles.	
ALMACENAMIENTO Y TRANSPORTE	Almacenar en lugar fresco y seco. Mantener alejado de la luz solar. Proteger de la humedad y evitar el calor excesivo.	
MSDS	No aplicable	

TECHNICAL DATA SHEET

BESTDENT MICRO APPLICATORS		
REFERENCES	25650: Thin 25651: Regular	
MANUFACTURER	CHL MEDICAL SOLUTIONS, Srl Via Felice Cavallotti, 15 20122 Milano - Italy	
CLASSIFICATION	Medical Device, Class I, according to Directive 93/42/EEC	
DESCRIPTION	The disposable micro applicator is an invasive device for transient use. It is used for the application of materials in the oral cavity. Sizes: Thin (yellow) and regular (green).	
COMPOSITION	Constituted by fiber micro brushes and a polypropylene stick.	
PACKAGING	Package with 100 pieces.	
SHELF LIFE	3 years from the date of manufacture	
INTENDED USE	The micro applicators are used for the application of small amounts of materials (etchants, bonding agents, sealers, cements) in limited access areas.	
SAFETY INFORMATION	For single use only. For professional dental application only. Non sterile.	
HANDLING AND STORAGE CONDITIONS	Store in a cool dry place, away from direct sunlight. Protect from moisture, and avoid excessive heat.	
MSDS	Not Applicable	