



FICHA TÉCNICA

PUNTAS DISPENSADORAS LUER LOCK BESTDENT		
REFERENCIAS	60019	
FABRICANTE / RESPONSABLE	CHL MEDICAL SOLUTIONS, Srl Via Felice Cavallotti, 15 20122 Milán - Italia	
CLASIFICACIÓN	Producto Sanitario Clase I, Directiva 93/42/CEE	
DESCRIPCIÓN	Las Puntas Dispensadoras son productos de uso en odontología para la limpieza de boca. Son curvadas 30° para facilitar el acceso a sitios complejos.	
COMPOSICIÓN	Las puntas están fabricadas de polipropileno PP sin ftalatos.	
ENVASE	Envases conteniendo 12 unidades	
VIDA ÚTIL	Debido a las características del producto, no tienen fecha prevista de caducidad.	
FINALIDAD PREVISTA	Las Puntas Dispensadoras son productos de uso en odontología para la aplicación de geles o líquidos en los dientes para la limpieza de la boca.	
INFORMACIÓN SOBRE LA SEGURIDAD	Solamente para uso profesional dental. Producto no estéril. No se considera un producto peligroso para la salud humana ni para el medio.	
ALMACENAMIENTO Y TRANSPORTE	Almacenar en lugar fresco y seco. Mantener alejado de la luz solar. Proteger de la humedad y evitar el calor excesivo.	

TECHNICAL DATA SHEET

BESTDENT LUER LOCK DISPENSING TIPS		
REFERENCE	60019	
MANUFACTURER	CHL MEDICAL SOLUTIONS, Srl Via Felice Cavallotti, 15 20122 Milano - Italy	
CLASSIFICATION	Medical Device, Class I, according to Directive 93/42/EEC.	
DESCRIPTION	The dispensing tips are products for professional dental use only used during the process of dental cleaning. They are bended 30° for easier access to the teeth.	
COMPOSITION	The tips are made of phthalate-free polypropylene.	
PACKAGING	Packages containing 12 pieces	
SHELF LIFE	Due to the characteristics of the product, no specific shelf-life is defined	
INTENDED USE	The dispensing tips are medical devices for the application of gels or liquids during dental cleaning.	
SAFETY INFORMATION	For professional dental use only. Non-sterile product. It is not considered a dangerous product for health or the environment.	
HANDLING AND STORAGE CONDITIONS	Store in a cool dry place, away from direct sunlight. Protect from moisture, and avoid excessive heat.	