

FICHA TÉCNICA

HILO DE RETRACCIÓN PROCLINIC		
REFERENCIAS	78287: Nº 000 78288: Nº 00 78289: Nº 0 78290: Nº 1	
FABRICANTE / RESPONSABLE	CHL MEDICAL SOLUTIONS, Srl Via Felice Cavallotti, 15 20122 Milán - Italia	
CLASIFICACIÓN	Producto Sanitario Clase I, según la Directiva 93/42/CEE	
DESCRIPCIÓN	El Hilo de Retracción es un hilo no impregnado trenzado para evitar que despeluchen, en los grosores más comunes del mercado. Son de fácil colocación y rápida retracción.	
COMPOSICIÓN	Fabricado en algodón 100 %.	
ESPECIFICACIONES	Algodón puro absorbente de dos líneas. Tamaños y colores: • 000: 0.5mm ± 0.05. Color negro. • 00: 0.6mm ± 0.05. Color amarillo. • 0: 0.8mm ± 0.05. Color morado. • 1: 0.9mm ± 0.05. Color azul.	
ENVASE	Envases con 2 metros de hilo.	
VIDA ÚTIL	Debido a las características del producto, el hilo de retracción no tiene una fecha de caducidad definida.	
FINALIDAD PREVISTA	Es un producto de uso en odontología para la retracción de encías y se suelen usar combinados con líquidos de retracción	
INFORMACIÓN SOBRE LA SEGURIDAD	Solamente para uso profesional dental. Producto de un solo uso. No estéril.	
ALMACENAMIENTO Y TRANSPORTE	Almacenar en lugar fresco y seco. Mantener alejado de la luz solar. Proteger de la humedad y evitar el calor excesivo.	

TECHNICAL DATA SHEET

RETRACTION CORD PROCLINIC		
REFERENCES	78287: No. 000 78288: No. 00 78289: No. 0 78290: No. 1	
MANUFACTURER	CHL MEDICAL SOLUTIONS, Srl Via Felice Cavallotti, 15 20122 Milano - Italy	
CLASSIFICATION	Medical Device, Class I, according to Directive 93/42/EEC	
DESCRIPTION	Non-impregnated braided retraction cord made of 100% cotton. Available in the most common sizes. Easy location and removal.	
COMPOSITION	Made of 100% cotton.	
SPECIFICATIONS	Absorbent cotton in two lines. Sizes and colours: • 000: 0.5mm ± 0.05. Black. • 00: 0.6mm ± 0.05. Yellow. • 0: 0.8mm ± 0.05. Violet. • 1: 0.9mm ± 0.05. Blue.	
PACKAGING	Package containing 2 meters of cord.	
SHELF LIFE	Due to the characteristics of the product, an expiration date is not defined.	
INTENDED USE	Retraction cord is used for rapid tissue displacement prior to impression making and can also be used together with retraction liquids.	
SAFETY INFORMATION	For single use only. For professional dental application only. Non sterile product.	
HANDLING AND STORAGE	Store in a cool dry place, away from direct sunlight. Protect from moisture, and avoid excessive heat.	