

	TECHNICAL DATA SHEET	CODE: TDS_00X
	DENTAL ROOT CANAL INSTRUMENTS	DATE: 09/05/18
	BESTDENT	VERSION: 01

PRODUCT	See annex I	
REFERENCE	See annex I	
COMPANY	SHENZHEN PERFECT MEDICAL INSTRUMENTS CO., LTD Add : 1F, Bldg. C Plant, No.3, Jixia Zaohekeng Industrial Area, Nanwan St., Longgang Dist., Shenzhen, China	
CLASSIFICATION	CLASS IIa MEDICAL DEVICE Directive 93/42/CEE	
PRODUCT DESCRIPTION	DENTAL ROOT CANAL INSTRUMENTS	
SPECIFICATIONS	File: stainless steel 304H, NITI. Handle material is divided into two kinds, one kind is plastic for hand use, another kind is stainless steel for engine use. Stopper: material is made of Silicone, and the material is Medical level.	
PACKAGING	Product is packed with PP plastic, Ten boxes into a middle carton. Carton: 5*33mm;10*33mm;5*28mm; Material: 275 lb. test double wall corrugated cardboard	
SHELF LIFE	Validity is defined as: 5 years. The Root canal files have been proven to be secure for tens of years.	
INTENDED OF USE	Root canal files are used for diagnostics, shaping/cleaning, filling in the root canal cure and for preparation of root canal filling.	
SECURITY	Root Canal instrument is considered to be used by a professional dentist in a dental clinic, in the patient's root canal treatment process according to the practical requirement.	
TRANSPORT AND STORAGE CONDITIONS	The metal material is very stable and does not be revised under conventional storage conditions.	

N/A: Not Aplicable

	TECHNICAL DATA SHEET	CODE: TDS_00X
	DENTAL ROOT CANAL INSTRUMENTS	DATE: 09/05/18
	BESTDENT	VERSION: 01

ANNEX I

87895 LIMAS PLURI TAPER STARTERKIT 21mm 6u BD
 87896 LIMAS PLURI TAPER STARTERKIT 25mm 6u BD
 87897 LIMAS PLURI TAPER STARTERKIT 31mm 6u BD
 87898 LIMAS PLURI TAPER SX 19mm 6u.BESTDENT
 87899 LIMAS PLURI TAPER S1 21mm 6u.BESTDENT
 87900 LIMAS PLURI TAPER S1 25mm 6u.BESTDENT
 87901 LIMAS PLURI TAPER S1 31mm 6u.BESTDENT
 87902 LIMAS PLURI TAPER S2 21mm 6u.BESTDENT
 87903 LIMAS PLURI TAPER S2 25mm 6u.BESTDENT
 87904 LIMAS PLURI TAPER S2 31mm 6u.BESTDENT
 87905 LIMAS PLURI TAPER F1 21mm 6u.BESTDENT
 87906 LIMAS PLURI TAPER F1 25mm 6u.BESTDENT
 87907 LIMAS PLURI TAPER F1 31mm 6u.BESTDENT
 87908 LIMAS PLURI TAPER F2 21mm 6u.BESTDENT
 87909 LIMAS PLURI TAPER F2 25mm 6u.BESTDENT
 87910 LIMAS PLURI TAPER F2 31mm 6u.BESTDENT
 87911 LIMAS PLURI TAPER F3 21mm 6u.BESTDENT
 87912 LIMAS PLURI TAPER F3 25mm 6u.BESTDENT
 87913 LIMAS PLURI TAPER F3 31mm 6u.BESTDENT
 87914 LIMAS PLURI TAPER F4 21mm 6u.BESTDENT
 87915 LIMAS PLURI TAPER F4 25mm 6u.BESTDENT
 87916 LIMAS PLURI TAPER F4 31mm 6u.BESTDENT
 87917 LIMAS PLURI TAPER F5 21mm 6u.BESTDENT
 87918 LIMAS PLURI TAPER F5 25mm 6u.BESTDENT
 87919 LIMAS PLURI TAPER F5 31mm 6u.BESTDENT
 87920 PLURI TAPER RETRAT. D1 (16mm.) 6u. BD
 87921 PLURI TAPER RETRAT. D2 (18mm.) 6u. BD
 87922 PLURI TAPER RETRAT. D3 (22mm.) 6u. BD
 87923 PLURI TAPER RETRAT. SURTIDO 6u. BD
 87924 LIMAS REC R25 - 21mm (6u.) BD
 87925 LIMAS REC R40 - 21mm (6u.) BD
 87926 LIMAS REC R50 - 21mm (6u.) BD
 87927 LIMAS REC SURTIDO -21mm (6u.)BD
 87928 LIMAS REC R25 - 25mm (6u.) BD
 87929 LIMAS REC R40 - 25mm (6u.) BD
 87930 LIMAS REC R50 - 25mm (6u.) BD
 87931 LIMAS REC SURTIDO -25mm (6u.) BD
 87932 LIMAS REC R25 - 31mm (6u.) BD
 87933 LIMAS REC R40 - 31mm (6u.) BD
 87934 LIMAS REC R50 - 31mm (6u.) BD
 87935 LIMAS REC SURTIDO -31mm (6u.) BD