


CHL MEDICAL SOLUTIONS, SRL

TECHNICAL DATA SHEET	Code: TDS_022
BESTDENT ARTICULATING PAPER	Date: 12/Ene/17
CHL	Version: 01

PRODUCT	ARTICULATING PAPER BESTDENT	
REFERENCE	01951 01956	
COMPANY	CHL MEDICAL SOLUTIONS SRL VIA FELICE CAVALLOTTI, 20122 MILANO - ITALY	
CLASSIFICATION	CLASS I MEDICAL DEVICE Directive 93/42/CEE	
PRODUCT DESCRIPTION	Paper with different thickness, forms, double-sided one color, and double-sided two colour.	
SPECIFICATIONS	<p>-Dimensions: 40µ or 200µ 22mmx110mm or 18mmx50mm</p> <p>-Appearance: Smooth cutting edges Nearly no protruding fibres Homogeneous intensity of colour Homogeneously coiled</p> <p>-Colour double-sided: red/red, blue/blue, red/blue</p>	
PACKAGING	Packed in a dispenser box made of Printed cardboard with front label (40µ) or white dispenser (plastic box) with front label (200µ) 120 sheets or 300 sheets: 100 per compartment	
SHELF LIFE	3 years	
INTENDED OF USE	Identification of uneven areas of dental restorations.	
SECURITY	This product is well known by the user and safe in use.	
TRANSPORT AND STORAGE CONDITIONS	Storage below 25 ° C	