



PRODUCT DATA SHEET

CHEMICAL/TECHNICAL NAME: Zinc Phosphate Cement

COMMERCIAL NAME: Improved Zinc Cement Type I
Powder 32 g + Liquid 15 mL

COMPOSITION:

- Zinc Oxide
- Magnesium Oxide
- Bismuth Oxide
- Silica Flour
- Fluoride Agent
- Hydrated Aluminum Phosphate
- Orthophosphorus
- Bidistilled Water

PHYSICAL FEATURES:

Pale, homogeneous yellow powder
Pale yellow liquid

USE AND FUNCTION:

- For cementing inlays, crowns and orthodontic bands.
- If the cavity is deep enough to compromise the pulp, it is advisable to use a layer of calcium hydroxide before using this cement.

DANGERS: May cause mild allergic reactions.