



DISPOSABLE DENTAL TRAY 28x18 COLORED - ROUGH BOTTOM

REFERENCE	BD4182.XX	CLIENTE:	
PACKAGING:	CONVENCIONAL	Customer ref:	
DRAWING:	<input checked="" type="checkbox"/> ADJUNTO		

IMAGE :



MATERIAL :

P.E.T.
THICKNESS : 180 +/- 5% MICRONS

COLOR :

DARK BLUE, CYAN, ORANGE, YELLOW

AVERAGE UNIT WEIGH:	14,50	grs.	TECHNIQUE: :	THERMOFORMED
---------------------	-------	------	--------------	--------------

PACKAGING:

EMBALAJE :	LOGO :	PROPIO CLIENTE	DISTRIBUTION :
	BOX MODEL:	14 click coin blanca	
	SIZE:	400X300X110 mm	
QUANTITY x BOX :	400 UNITS		2 piles of 200 units each, in a unique polyethylene bag

PALETIZATION:

PALLET:	EUROPALET 800 X 1200
MAXIMUM HEIGH:	1.780,00 mm
DISTRIBUTION:	18 boxes x 5 stages + 6 boxes
QUANTITY x PALLET:	38.400 UNITS
FINAL WRAPPING:	Shrink-Wrapping
Packaging observations:	

REGULATION:

Medical device class 1, European directive 93/42 CEE

QUALITY DEPT.	CLICK COIN, S.L.	CUSTOMER COMPLIANCE	
QUALITY MANAGER	Gemma Bochaca Royo	RESPONSIBLE:	
DATE:	14-05-18	DATE:	
SIGNATURE:		SIGNATURE:	(If needed)