


## 1. Products Description

### 1.1 Products Name: Bouffant cap

### 1.2 Products Details:

Dochem Code	Description	Package	Products Photo
1A2700	18''(46cm), white	100 pcs/bag, 20 bags/ctn	
1A2701	18''(46cm), green	100 pcs/bag, 20 bags/ctn	
1A2702	18''(46cm), blue	100 pcs/bag, 20 bags/ctn	
1A2703	18''(46cm), pink	100 pcs/bag, 20 bags/ctn	
1A2704	21''(51cm), white	100 pcs/bag, 20 bags/ctn	
1A2705	21''(51cm), green	100 pcs/bag, 20 bags/ctn	
1A2706	21''(51cm), blue	100 pcs/bag, 20 bags/ctn	
1A2707	21''(51cm), pink	100 pcs/bag, 20 bags/ctn	

### 1.3 Materials of Products

Bouffant caps are made of PP Spunbonded Non-woven Fabrics (Polypropylene)

which are ruled by relevant standards of Europe.

### 1.4 Purpose of Products

Bouffant caps are used to help dentists to dust proof, anti-hair loss during treatment to keep a clean treatment environment.

### 1.5 Usages of Products

1.5.1 Only dentists and nurse can use it.

1.5.2 Choose the best right size according to the user's size.

### 1.6 Notice of Use

1.6.1 Bouffant caps are one time use, can not be repetition used.

1.6.2 If the product is broke, please don't use.

## 1.7 Storage Condition

1.7.1 Keep in the well ventilated rooms, free from corrosive air

1.7.2 Temperature control  $\leq 40C$ , moisture  $\leq 80\%$ .

1.8 Years of Shelf Life: 3 years

## 1.9 International Products Code

According to section I definition and section III assortment of MDD93/42EEC:

its classification: **I**

Product Name	Time of Use	Plug-in Medical Device	No. of Usage	Sterilization	Classification	UMDNS Code
Bouffant caps	Less than 60 minutes	no	Single Use	no	I	16081


## 2. Sterilization Methods

This product do not need sterilize.

## 3. Labels and Language

3.1 Company Information: Details of the manufacturer and EU representative.

### 3.1.1 Manufacture

a) Manufacture symbol: 

b) Manufacture information:

Shanghai Dochem Industries Co., Ltd.

Address: Suite E2014, 2nd Floor, No. 3558, ZhenBei Rd. Putuo district, Shanghai 200062, China.

info@dochem.biz dochem@sh163.com

www.dochem.biz

### 3.1.2 EU representative

a) EU representative: 

EC	REP
----	-----

b) EU representative information:

MedNet EC-Rep GmbH

Borkstrasse 10, 48163 Muenster, Germany

### 3.2 Products Information: Product name, type, quantity and specification.

### 3.3 All Signs:

3.3.1 Sterilization sign:



3.3.2 Lot System:

a) Lot symbol: 

LOT
-----

B) Lot number: **XX XX XXXX**

**(XX: year of produce, XX: month of produce, XXXX: our contract number)**

Example: 18 0 001 (Production year 2018, no meaning 0, contract number 001)

3.3.3 Manufacture Date:


a) Manufacture symbol: 

b) Manufacture Number: **XXXX XX XX**

**(XXXX: year of produce, XX: month of produce, XX: date of produce)**

Example: 20180101 (Produce year 2018, month 01, date 01)


**3.3.4 Expiration Date:**

a) Expiration symbol: 

b) Expiration Number: **XXXX XX XX**

**(XXXX: year of expiry, XX: month of expiry, XX: date of expiry)**

Example: 20170101 (Produce year 2017, month 01, date 01)

3.3.5 One time use symbol: 

3.3.6 Guidance instruction symbol: 

3.3.7 Package broken symbol: 

3.3.8 Made in China