

TECHNICAL DATA SHEET

**SPATULA HEIDEMANN "1045/12B"
PROCLINIC**



REFERENCES	flexible - 2,5mm
MANUFACTURER	Carl Martin GmbH Neuenkamper Str. 80-86 42657 Solingen, Germany
CLASSIFICATION	Medical Devices class I, according to Directive 93/42/EWG
DESCRIPTION	Spatula Heidemann
COMPOSITION	surgical stainless steel
SPECIFICATIONS	Autoclavable up to 134°C. Resistant to the chemical sterilant. Supplied in "Non-sterile" condition
PACKAGING	1 unit
INTENDED USE	This instrument is used in conservative restorative therapy. It is used for the introduction, plugging, compaction, smoothing and modeling of filling materials (cement, composite, amalgam etc.)
SHELF LIFE	Not applicable
SAFTEY INFORMATION	Intended for use by dental professionals only. Must be checked for fault-free operation (fractures, cracks, malfunctions, etc.) and be cleaned, disinfected and sterilized before each use. New instruments are contaminated with oils and fats from the manufacturing process and must therefore be cleaned thoroughly before the first sterilisation with a suitable cleaning agent.
HANDLING AND STORAGE	Sterile products must be stored in a dry, clean and dust-free environment at temperatures between 5°C and 40°C.

Solingen, 01.10.2018