



PRODUCT DATA SHEET

CHEMICAL/TECHNICAL NAME: Cement Z.O.E Type III

COMMERCIAL NAME: Cement Z.O.E Type III Intermediate Restorative Material
Powder 38 g + Liquid 14 mL

COMPOSITION:

- Zinc Oxide
- Aluminum Oxide
- Acrylic Polymer
- Eugenol
- Hardener

PHYSICAL FEATURES:

Pale, homogeneous white powder
Pale yellow to amber liquid

USE AND FUNCTION:

- Intermediate restorations particularly in pediatric orthodontics and decayed teeth.
- Protection of the tooth after removing demineralized tissues from the cavities.
- It can also be used for bases or as temporary cement, it also has the advantage of protecting and healing the pulp.

DANGERS: May cause mild allergic reactions.