

## MOUTH MIRRORS DOUBLE SIDED

| Steel           | Aspect | Origin |
|-----------------|--------|--------|
| Stainless Steel | Shiny  | France |

| Reference | Medical Device Class 1          |
|-----------|---------------------------------|
| MPP8324PH | PLANE N 5 DOUBLE FACE CSS BOX 6 |

| stem steel               | simple stem size                             | Frame steel         |
|--------------------------|----------------------------------------------|---------------------|
| simple stem: 303 EN10088 | Ø2,5mm - Length 32mm - Fillet ISO 2,5x0,45mm | 316L EN10088:1.4404 |

| stem steel                            | cone socket stem size                       | Silicone glue |
|---------------------------------------|---------------------------------------------|---------------|
| Cone socket stem: 316L EN10088:1.4404 | Ø2,84mm - Length 50mm - Fillet USA 2,51x56" | RTV732 Clear  |

| Mirror glass                                                     | Size                                                                         |
|------------------------------------------------------------------|------------------------------------------------------------------------------|
| Aluminium/titanium/Pure reflect through evaporation under vacuum | Thickness 1,6/1,8mm - Ø22/24mm according to mirror size - Tolerance +-0,10mm |

### Marking

#### Identification

Manufacturer *Prodont-Holliger*  
Item code  
Lot number

**Marking:** Laser  
CE

#### Cleaning

All means  
*except metallic or abrasive*  
(Ex. Scotch brite)

#### Sterilisation

all (autoclave, ...)  
*except corrosive products*  
follow cleaning products manufacturers instructions

#### Warranty

*exchange if defect, except normal wear*

#### Storage

Dry, clean place

#### Shipping

48 hours on stock

#### Maintenance

N/A

#### Medical product surveillance resource

Audrey Vialle

Last revised version: 2018/01